# Operating Instructions

CO<sub>2</sub> Incubator

**MCO-170AIC**

**MCO-170AICUV**

**MCO-170AICUVH**

Please read the operating instructions carefully before using this product, and keep the operating instructions for future use.

See page 72 for model numbers.

# CONTENTS

INTRODUCTION .....	4
INTENDED USE AND PRECAUTIONS.....	4
PRECAUTIONS FOR SAFE OPERATION.....	5
LABELS ON INCUBATOR .....	10
SYMBOLS ON INCUBATOR .....	11
ENVIRONMENTAL CONDITIONS.....	11
INCUBATOR COMPONENTS .....	12
Unit.....	12
LCD touch panel.....	14
Remote alarm terminal.....	16
INSTALLATION.....	17
Installation site.....	17
Installation .....	18
Connecting a CO <sub>2</sub> gas cylinder .....	21
BEFORE COMMENCING OPERATION.....	22
Initial cleaning method.....	22
Removing inner attachments .....	23
Installing inner attachments.....	25
Filling the humidifying pan.....	26
FOR OPTIMAL CULTIVATION.....	27
Precautions for cultures.....	27
Preventing contamination.....	28
CORRECT OPERATION .....	29
BASIC OPERATION ON LCD TOUCH PANEL .....	30
Operation from Menu key.....	30
Entering numeric values into input window .....	32
Operating slide keys.....	33
BASIC PARAMETERS.....	34
Setting temperature, CO <sub>2</sub> density, and high-limit temperature alarm .....	34
Setting key lock .....	35
Turning off key lock .....	35
ALARM PARAMETERS .....	36
OPERATION/ALARM LOG .....	38
Setting log interval.....	38
Displaying operation log.....	39
Exporting operation log .....	41
Displaying alarm log.....	43
Exporting alarm log .....	44
OTHER PARAMETERS.....	46
Setting date and time .....	46
Setting brightness and sleep.....	47

# CONTENTS

UV LAMP PARAMETERS.....	48
Using the UV lamp .....	48
Setting UV lamp ON period.....	49
Lighting the UV lamp for 24 hours.....	50
H <sub>2</sub> O <sub>2</sub> DECONTAMINATION .....	51
H <sub>2</sub> O <sub>2</sub> decontamination .....	51
Precautions when handling H <sub>2</sub> O <sub>2</sub> reagent.....	55
ELECTRIC LOCK (OPTIONAL).....	56
Setting User-ID .....	56
Setting auto-lock.....	57
Using unlock key .....	59
Turning off auto-lock.....	59
GAS AUTO-CHANGER (OPTIONAL).....	60
Connecting CO <sub>2</sub> gas cylinders .....	60
Automatic CO <sub>2</sub> gas supply line changeover .....	60
Manual CO <sub>2</sub> gas supply line A/B changeover .....	62
STD GAS AUTO-CALIBRATION KIT (OPTIONAL) .....	63
ALARMS, SAFETY, AND SELF-DIAGNOSIS.....	65
TROUBLESHOOTING.....	68
SPECIFICATIONS .....	70
PERFORMANCE .....	72
SAFETY CHECK SHEET .....	73

---

# INTRODUCTION

---

- Read the operating instructions carefully before using the product and follow the instructions for safe operation.
- PHC Corporation takes no responsibility for safety if the product is not used as intended or is used with any procedures other than those given in the operating instructions.
- Keep the operating instructions in a suitable place so that they can be referred to as necessary.
- The operating instructions are subject to change without notice for improvement of performance or function.
- Contact our sales representative or agent if any page of the operating instructions is lost or the page order is incorrect, or if the instructions are unclear or inaccurate.
- No part of the operating instructions may be reproduced in any form without the express written permission of PHC Corporation.

## IMPORTANT NOTICE

PHC Corporation guarantees this product under certain warranty conditions. However, please note that PHC Corporation shall not be responsible for any loss or damage to the contents of the product.

---

# INTENDED USE AND PRECAUTIONS

---

This equipment is designed for culture of cells, tissues, organs, and embryos.

- The adapted culture condition depends on the kind of sample. It is necessary to determine the culture temperature, CO<sub>2</sub> density, and culture period that is required.
- For culture of embryos, low O<sub>2</sub> density is desirable. Using an O<sub>2</sub>/CO<sub>2</sub> incubator is recommended.
- For IVF/ART use, special attention should be paid to traceability, since incidents may be disclosed several months or years later, at the baby's birth, or during its later life. It is therefore recommended that the following data is maintained: product serial number, incubation term, and incubation parameters. (Refer to details on MEDDEV 2.2/4.)

# PRECAUTIONS FOR SAFE OPERATION

**It is imperative that the user complies with the operating instructions as they contain important safety advice.**

Items and procedures are described so that this unit can be used correctly and safely. Following these precautions will prevent possible injury to the user and any other person.

Precautions are illustrated in the following way:

## **WARNING**

Warning indicates a potentially hazardous situation which, if not avoided, could result in serious injury or death.

## **CAUTION**

Failure to observe CAUTION signs could result in injury to personnel and damage to the unit and associated property.

Symbols have the following meaning;

-  This symbol means caution.
-  This symbol means an action is prohibited.
-  This symbol means an instruction must be followed.

Be sure to keep the operating instructions in a place that is accessible to users of this unit.

## **WARNING**

As with any equipment that uses CO<sub>2</sub> gas, there is a likelihood of oxygen depletion in the vicinity of the equipment. It is important that you assess the work site to ensure there is suitable and sufficient ventilation. If restricted ventilation is suspected, then other methods of ensuring a safe environment must be considered. These may include atmosphere monitoring and warning devices.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  **Do not use the unit outdoors.** Exposure to rain may cause leakage and/or electric shock.
-  **Only qualified engineers or service personnel should install the unit.** The installation by unqualified personnel may cause electric shock or fire.
-  **Install the unit in a location capable of bearing the total combined weight (product + optional accessories + stored items). After installing the unit, be absolutely sure to take precautions to prevent the unit from falling over.** If the unit is installed in a location which is not strong enough or if the proper precautions are not taken, the unit may fall over and cause injuries.
-  **Do not install the unit where there are high levels of moisture or where it may be splashed with water.** Installing the unit where there are high levels of moisture or where it may be splashed with water may cause the insulation to deteriorate and give rise to leakage and/or electric shock.
-  **Do not install the unit in a location where flammable or volatile substances are present.** Installing the unit in a location where flammable or volatile substances are present may cause explosions and/or a fire.
-  **Do not install the unit in a location where corrosive gases such as acids are present.** Installing the unit in a location where corrosive substances are present may cause electric components to corrode, leading to leakage and/or electric shock due to the deterioration of insulation resulting from corroded electrical components.
-  **Do not place this unit in a location where it is difficult to disconnect the power supply plug.** Failure to disconnect the power supply plug may cause fire in the event of a problem or malfunction.
-  **Be absolutely sure to earth (ground) the unit in order to prevent electric shock.** Failure to earth the product may give rise to electric shock. If necessary, ask a qualified contractor to do this work.
-  **Do not connect the earth wire to a gas pipe, water pipe, or lightning rod when earthing the unit.** Earthing the unit improperly may give rise to electric shock.
-  **Connect the unit to a power source as indicated on the rating label attached to the unit.** Use of any other voltage or frequency other than that on the rating label may cause fire or electric shock.
-  **Never store volatile or flammable substances in this unit except in a sealed container.** Such substances may cause explosion or fire if they leak.
-  **Never insert metal objects such as pins and wires into any vent, gap or outlet on the unit.** This may cause electric shock or injury by accidental contact with moving parts.
-  **When handling harmful samples (for example, those which consist of toxic, pathogenic, or radioactive substances), install the unit inside a designated isolation facility.** If the unit is installed in a location which is not an isolation facility, there may be detrimental effects on both people and the natural environment.

# WARNING

-  **Before proceeding with maintenance or checking of the unit, set the power switch to OFF and disconnect the power supply plug.** Performing the work while power is still flowing to the product or while the power-supply plug is still connected may give rise to electric shock and/or injury.
-  **Do not touch any electrical parts (such as power supply plug) or operate switches with a wet hand.** This may cause electric shock.
-  **Wear protective gloves and mask during maintenance.** Touching or inhaling chemicals or aerosols from around the unit may be detrimental to health.
-  **Never splash water directly onto the unit** as this may cause electric shock or short circuit.
-  **Never put containers with liquid on top of the unit** as this may cause electric shock or short circuit if the liquid is spilled.
-  **Never damage the power supply cord or power supply plug (by breaking, adapting, placing near a source of heat, bending with force, twisting, pulling, adding weight, or binding).** A damaged power supply cord or power supply plug may cause electric shock, short circuit, or fire.
-  **Never disassemble, repair, or modify the unit yourself.** A high-voltage area is located inside the unit. Any work carried out by an unauthorized personnel may result in electric shock. Contact our sales representative or agent for maintenance or repair.
-  **Make sure the power supply plug is pushed fully in.** Faulty insertion of the power supply plug may cause electric shock or fire due to generation of heat. Never use a damaged power supply plug or loose power outlet.
-  **Disconnect the power supply plug if there is anything wrong with the unit.** Continued abnormal operation may cause electric shock or fire.
-  **Grip the power supply plug when disconnecting the power supply cord from the outlet.** Pulling the power supply cord may cause electric shock or short circuit.
-  **Remove dust from the power supply plug periodically.** Dust on the power supply plug may cause insulation failure due to moisture and thus cause a fire. Disconnect the power supply plug and wipe it with a dry cloth.
-  **Disconnect the power supply plug** before moving the unit. Take care not to damage the power supply cord. A damaged cord may cause electric shock or fire.
-  **Disconnect the power supply cord when the unit is not in use for long periods.** Keeping the unit connected may cause electric shock, leakage, or fire due to the deterioration of insulation.

# PRECAUTIONS FOR SAFE OPERATION

## **WARNING**

-  If the unit is to be stored unused in an unsupervised area for a long period, **ensure that children do not have access and that doors cannot be closed completely.**
-  **Ask a qualified contractor to carry out disassembly and disposal of the unit.** Leaving the unit in a location that can be accessed by third parties may result in unexpected accidents (e.g. the unit may be used for unintended purposes).
-  **Do not leave the plastic bags used for packing in a place where they can be reached by small children** as this may result in unexpected accidents such as suffocation.
-  **Use the reagent specified by our company for H<sub>2</sub>O<sub>2</sub> decontamination.** Using a different H<sub>2</sub>O<sub>2</sub> solution may result in explosion or damage to the incubator.
-  When performing H<sub>2</sub>O<sub>2</sub> decontamination, **securely close the internal and external doors.** Failure to do so may be harmful to health due to leakage of H<sub>2</sub>O<sub>2</sub> gas.
-  During H<sub>2</sub>O<sub>2</sub> decontamination, **plug the access hole with the silicon cap that is provided.** Failure to do so may be harmful to health due to leakage of H<sub>2</sub>O<sub>2</sub> gas.
-  **Always use the removal power supply cord that is provided.** Other power supply cord may cause electric shock or fire.
-  When using CO<sub>2</sub> gas for control, **make sure that there is adequate ventilation.** Using CO<sub>2</sub> gas in a small room without adequate ventilation may cause gas poisoning or oxygen deprivation. In addition, when opening the incubator doors, do not directly inhale the air in the chamber.
-  When connecting a gas cylinder to the incubator, **confirm the gas type. Confirm that the connections are secure and that no gas will leak. Be sure to use the specified pressure.** Using an incorrect gas or pressure may result in explosion or fire, or in gas poisoning or oxygen deprivation due to gas leak.
-  **Install the incubator in a location with adequate ventilation.** If adequate ventilation cannot be provided, then install an alarm system using CO<sub>2</sub> and O<sub>2</sub> densitometers.
-  **Do not look directly at UV light.** UV light is harmful to the eyes.
-  Never start H<sub>2</sub>O<sub>2</sub> decontamination when “Warning: UV Bulb Life” is displayed in the message display field. The UV resolve is not sufficient.
-  **Do not use the unlock key to unlock the outer door during H<sub>2</sub>O<sub>2</sub> decontamination or during H<sub>2</sub>O<sub>2</sub> gas resolve by UV.** Doing so may cause harm to health from H<sub>2</sub>O<sub>2</sub> gas leakage.
-  **The CO<sub>2</sub> incubator must be dismantled and disposed of by qualified personnel only.** If the CO<sub>2</sub> incubator is left where outsiders enter, it may result unexpected accident (for example, children to become locked inside).

# CAUTION

-  This unit must be plugged into a dedicated circuit protected by branch circuit breaker.
-  Use a dedicated power source as indicated on the rating label attached to the unit. A multiple-tap may cause fire resulting from abnormal heating.
-  **Never store corrosive substances such as acid or alkali** in this unit if the container cannot be sealed. These may cause corrosion of inner components or electric parts.
-  **Check the setting when starting up of operation after power failure or turning off of power switch.** The stored items may be damaged due to the change of setting.
-  **Be careful not to tip over the unit** during movement to prevent damage or injury.
-  **Prepare a safety check sheet** (copy the last page) when you request any repair or maintenance for the safety of service personnel.
-  **Wear rubber gloves when handling the H<sub>2</sub>O<sub>2</sub> reagent.** Direct contact with the H<sub>2</sub>O<sub>2</sub> reagent may result in inflammation of the skin.
-  H<sub>2</sub>O<sub>2</sub> decontamination can be performed only for the chamber and chamber attachments with standard specifications, and not for any other objects.
-  **Perform H<sub>2</sub>O<sub>2</sub> decontamination with the chamber attachments arranged as specified by our company.** Arranging them in a different way may result in insufficient decontamination.
-  After H<sub>2</sub>O<sub>2</sub> decontamination has been completed, **wear rubber gloves and use a non-woven cloth to wipe off the residual H<sub>2</sub>O<sub>2</sub> fluid from the bottom of the chamber, any objects that were decontaminated, and the bottoms of ducts.**
-  **Do not climb on top of the incubator or place any object on it.** By doing so you may damage the incubator or knock it over, which may result in injury.
-  **Wear rubber gloves when performing maintenance on the chamber.** Failure to wear gloves may result in cuts or abrasions from sharp edges or corners.

# LABELS ON INCUBATOR

Users are advised to read carefully the warnings and cautions contained on stickers at key locations on the interior and exterior of the incubator.

Possible Danger	Warning/Caution Type Location of Danger	Warning/Caution Label	Description of Danger
Burns	<b>Hot Surface</b> Cooling Unit & Heat Cover		Avoid touching the cooling unit and heat cover, which can reach high temperatures and may cause burns.
Personal injury	<b>Hazardous UV Light</b> Interior		The UV lamp is illuminated by pressing the door switch. Do not press the door switch: the UV light is hazardous.
Personal injury	<b>Hazardous UV Light</b> Interior		The UV light is hazardous. Never turn on the UV lamp without the cover.
Personal injury	<b>Gas Poisoning or Oxygen Deprivation</b> Environment		When using CO <sub>2</sub> gas, <b>ensure there is adequate ventilation</b> . Using CO <sub>2</sub> gas in a small room without adequate ventilation may cause gas poisoning or oxygen deprivation. In addition, when opening the incubator doors, do not directly inhale the air in the chamber.
Personal injury	<b>Gas Poisoning or Oxygen Deprivation</b> Interior		Excessive pressure may cause gas supply lines inside the incubator to come loose, which may result in gas poisoning or oxygen deprivation due to the escape of gas.
Personal injury	<b>Electric shock</b> Top cover		Never remove this cover. Only a service personnel remove the cover to avoid the electric shock.

# SYMBOLS ON INCUBATOR

The following symbols are attached to the incubator:

	Attached to covers that access high-voltage electrical components to prevent electric shock. Only a qualified engineer or service personnel should be allowed to open these covers.
	Indicates an ultraviolet light (UV) caution.
	Indicates that caution is required. Refer to precautions for safe operation for details.
	Indicates a hot surface.
	Indicates an earthing.
	Indicates "ON" for a power switch.
○	Indicates "OFF" for a power switch.




## ENVIRONMENTAL CONDITIONS

This equipment is designed to be safe at least under the following conditions (based on the IEC 61010-1):


- Indoor use;
- Altitude up to 2000 m;
- Temperature 5 °C to 40 °C;
- Maximum relative humidity 80 % for temperature up to 31 °C decreasing linearly to 50 % relative humidity at 40 °C;
- Mains supply voltage fluctuations up to  $\pm 10$  % of the nominal voltage;
- Transient overvoltages up to the levels of OVERVOLTAGE CATEGORY II;
- Temporary OVERVOLTAGES occurring on the mains supply;
- Applicable pollution degree of the intended environment (POLLUTION DEGREE 2 in most cases);

# INCUBATOR COMPONENTS

## Unit

(When some inner attachments are removed)


### Lower right side

(When 13 is removed)

Sample air outlet cap

Removable power supply cord port

\*1: MCO-170AICUVH/170AICUV or when an optional UV system set MCO-170UVS is installed.

\*2: MCO-170AICUVH or when an optional electric lock MCO-170EL is installed.

\*3: Type FG-7P or FG-7PL.


- 1. Outer door:** The outer door is held to the frame with the magnetic seal. The door heater is installed in the door panel. The door opening is reversible. Contact our sales representative or agent to change the door hinge from left to right or vice versa.
- 2. Inner door:** The inner door is made of tempered glass. Do not subject the glass to excessive impact. Contact our sales representative or agent to change the inner door to optional small doors (MCO-170ID).
- 3. Tray catches:** Insert tray to fit the concave portion on the chamber.
- 4. Fan cover (above the duct):** Serves as the inlet for circulating air. It is removable.
- 5. Duct:** The duct for the path for circulating air. It is removable.
- 6. Fan (inside the duct):** The fan is made from polypropylene resin. It can be disinfected in an autoclave.
- 7. UV lamp\*<sup>1</sup>:** This UV lamp does not generate ozone. For replacement, contact our sales representative or agent.
- 8. Humidity control bar:** Automatically reduces condensation caused by the external environment when the door is frequently opened and closed. The humidity control bar has a bactericidal effect due to its plated surface. However, it is recommended that it is replaced every 5 to 6 years to maintain this effect (The duration of bactericidal effect differs depending on the environment).
- 9. Humidifying pan cover:** Prevents UV light entering the chamber.
- 10. Humidifying pan:** Fill with sterile distilled water.
- 11. Door switch:** Detects the door opening/closing and stops the fan and electromagnetic valve for CO<sub>2</sub> when the door is open. The UV lamp\*<sup>1</sup> is also deactivated by the door opening.
- 12. Key hole\*<sup>2</sup>:** This is unlocked with the unlock key when the outer door is locked by the electric lock.
- 13. Switch cover:** Prevents shutting down of the unit in case of accidental pressing of the main power switch.
- 14. Power supply cord cover plate:** Prevents the power supply cord becoming detached.
- 15. Connecting port A for CO<sub>2</sub> gas pipe:** Refer to page 21 for gas cylinder connection.  
**Note:** When the optional MCO-21GC gas auto-changer is installed, both ports A and B are available. Refer to page 60.
- 16. Glow starter\*<sup>3</sup>:** Starts the glow for the UV lamp.
- 17. Power switch:** This is the main switch for the incubator (ON-"I", OFF-"O"). It also functions as an overcurrent breaker.
- 18. Sample air outlet:** This also functions as an internal gas outlet. In normal use, cover this outlet with the sample air outlet cap.
- 19. USB port:** Insert USB memory device to export operations and alarms log. Refer to pages 41–45.
- 20. Access port:** Place the silicon caps on both the outside and inside of the port when the port is not being used.
- 21. Remote alarm terminal:** Sends the alarm to a remote location by connecting to an external alarm unit. Refer to page 16.

# INCUBATOR COMPONENTS

## LCD touch panel

The following display (called the Top screen) will appear when the power switch is turned ON.

**Note:** It takes approximately 20 seconds until the Top screen is displayed. During warming-up, "Status: Gas sensor initializing" is displayed in the Message display field (13), and "--." is displayed in the present CO<sub>2</sub> density display field (4).

### 1. Present temperature display field

The current chamber temperature is displayed.

### 2. Set temperature value display field

The set value of chamber temperature is displayed. Default setting: 37 °C

### 3. Heating indicator

This lights up when the heater is on.

### 4. Present CO<sub>2</sub> density display field

The current chamber CO<sub>2</sub> density is displayed. Nothing is displayed when CO<sub>2</sub> density is set to 0 %.

### 5. Set CO<sub>2</sub> density value display field

The set value of the chamber CO<sub>2</sub> density is displayed. Default setting: 0 %.

### 6. CO<sub>2</sub> gas injection indicator

This lights up when CO<sub>2</sub> gas is being injected.


### 7. CO<sub>2</sub> gas supply line indicator A and B\*

Current supplying CO<sub>2</sub> gas supply line (connecting port for CO<sub>2</sub> gas pipe) is displayed. When the CO<sub>2</sub> gas automatic changer function changes the empty CO<sub>2</sub> gas cylinder to the other, the empty indicator is displayed in reverse video and blinks.

### 8. CO<sub>2</sub> gas supply line select key\*

This key selects CO<sub>2</sub> gas supply line A or B (connecting port A or B for CO<sub>2</sub> gas pipe). When an optional gas auto-changer MCO-21GC is installed, CO<sub>2</sub> gas supply line A/B changes over automatically when CO<sub>2</sub> gas cylinder is empty.

\* Only displayed when optional component MCO-21GC (gas auto-changer) is installed.

### 9. Present date/time display field

The date and time is set when the incubator is shipped from the factory. Refer to page 46 for details.

### 10. Over heat display

High-limit temperature alarm is activated: "Over Heat" is displayed alternately in normal characters and reverse video.

### 11. UV lamp condition display

UV lamp ON: "UV : On" is displayed.  
UV lamp OFF: "UV : Off" is displayed.

**Note:** Displayed in models MCO-170AICUVH/170AICUV or when the optional UV system set MCO-170UVS is installed on the model MCO-170AIC.

### 12. Outer door (opening/closing) display

Open: "Door : Open" is displayed alternately in normal characters and reverse video.  
Closed: "Door : Closed" is displayed.  
Locked: "Door : Locked" is displayed.

### 13. Message display field

Alarms, errors, or messages are displayed when a fault occurs. Refer to pages 65–67.

**Note:** When there are several errors, the display shows the message number. For example, if there are two errors, the display shows "1/2".

### 14. Message select key

When there are several alarms/errors, the message on the screen can be changed.

### 15. Menu key

Press this key to access the Menu screen. Refer to page 30.

### 16. H<sub>2</sub>O<sub>2</sub> key

This key initiates H<sub>2</sub>O<sub>2</sub> decontamination. Refer to pages 51–55. If the incubator is not equipped for H<sub>2</sub>O<sub>2</sub> decontamination, this key is not displayed on the LCD touch panel.

### 17. Unlock key

Press this key to unlock the outer door when it is auto-locked by the electric lock. Refer to page 58. If the incubator is equipped for auto-lock and the auto-lock function is ON, this key is displayed.

### 18. Buzzer key

Press this key to silence the buzzer. Refer to pages 37 and 65–67.

# INCUBATOR COMPONENTS

## Remote alarm terminal


The alarm of this product can be transferred to a remote location by connecting the external alarm unit to the remote alarm terminal. For the behaviour of remote alarm output, refer to pages 65–67.

The terminal of the remote alarm is installed at the rear upper right of the unit (see the figure opposite). The alarm outputs from this terminal. Contact capacity is DC 30 V, 2 A.

Table 1 shows the behaviour of the remote alarm when the Buzzer key is pressed.

**Notes:**

- The door alarm does not work remotely. Refer to pages 65–66.
- Wiring of the remote alarm should be performed by qualified service personnel.

**Table 1. Behaviour of remote alarm when pressing Buzzer key**

Remote Alarm setting (Refer to pages 36–37)	Connecting terminal	Normal condition	Abnormal condition (Including in the cases of power outage and of where the power supply plug is pulled out)	
			When pressing the buzzer key	
ON: Remote alarm setting not connected with Buzzer key (Remote alarm is cancelled)	COM.-N.C.	Close	Open	Open (Maintain in abnormality)*
	COM.-N.O.	Open	Close	Close (Maintain in abnormality)*
OFF: Remote alarm setting connected with Buzzer key (Remote alarm is not cancelled)	COM.-N.C.	Close	Open	Close (Return to normal)
	COM.-N.O.	Open	Close	Open (Return to normal)

\* In the case of Err01 (CO<sub>2</sub> gas cylinder empty), Err11 and Err12 (CO<sub>2</sub> sensor error), the condition returns to normal.

Use a twisted shield wire for the connection.  
 Type: UL 2343, UL 2448, UL 2464, UL 2552, UL2623.  
 Length: 30 m max.

# INSTALLATION

## Installation site

For correct operation, install the incubator in a location with the following conditions (referring also to the precautions for safe operation on pages 5–9):

- **Install in an environment with normal air.**
- **Do not expose to direct sunlight.**
- **Do not install near significant heat sources, such as heaters, boilers, ovens, or autoclaves.**
- **Ensure ambient temperature is at least 5 °C lower than set temperature.**
- **Do not allow the ambient temperature to become too high.**
- **Install in a site with a relative humidity of 80 %R.H. or lower.**
- **Do not expose to flammable or corrosive gas.**
- **Do not install where objects may fall on the incubator from above.**

# INSTALLATION

## Installation


### 1. Remove the packing tape and clean up.

Remove all the tapes that are securing the doors and the inner attachments. Open the doors for ventilation. If the outer panels are dirty, use a cloth to wipe them with a diluted neutral detergent (Undiluted detergent can damage the plastic components. For the dilution, refer to the instructions on the detergent). Wipe off the residual detergent with a wet cloth and then wipe off any moisture.

**Note:** Remove the cable tie around the power supply cord to prevent corrosion of the cord coating.


### 2. Set up the humidifying pan and humidifying pan cover (Fig. 1).

### 3. Install 4 trays (Fig. 2).

(When optional small doors (MCO-170ID) are installed)


To use the gastight split doors effectively, set 3 trays at the same level as shown in Fig 3 below.

#### 4. Adjust the levelling feet.

Adjust the levelling feet by turning them anticlockwise until the incubator is level (Fig. 4).

**Note:** Incubating on a sloping tray may adversely affect the cultivation.

#### 5. Earth the incubator.

Earth the incubator during installation to prevent electric shock. If there is no earth wire at the location, consult with qualified service personnel.

- **When an earth must be installed**

If an earthed 3-pole outlet is not available, then an earth must be installed. Consult with qualified service personnel.

- **Installing an earth fault circuit breaker**

If you must use the incubator in a moist or humid location, then it is recommended that an earth fault circuit breaker be installed in the power supply circuit (i.e., the power supply at the incubator). Have the circuit breaker installed by qualified service personnel.


# INSTALLATION

## ● Double-stacking

For stacking the incubators securely, refer to the procedure included with the optional double-stacking bracket MCO-170PS or the stacking plate MCO-170SB.

**Note:** Two hooks are attached at the upper rear of the incubator. When stacking incubators, fix the upper incubator to the wall with these hooks and wire or chain (Fig. 5).

**Note:** When stacking the incubators with our CO<sub>2</sub> incubator or another O<sub>2</sub>/CO<sub>2</sub> incubator, use the stacking plate MCO-170SB. Refer to Table 11 on page 71.

## ● When the incubator is not in use

Empty the water from the humidifying pan and remove moisture from the chamber. Make sure that the chamber is completely dry before closing the doors. Failure to do so may result in damage.

## ● Before moving the incubator

Before moving the incubator, empty the water from the humidifying pan, disconnect the power supply plug from the outlet, and make sure that the power supply cord will not be damaged. Failure to do so may result in electric shock or fire.

## ● If the key is fitted in the door lock

Remove the key from the door lock and keep it in a secure place. Make a note of the key symbol and number in case the key is lost. To remove the key from the door lock, open the outer door and turn the key in the direction of the lock. Remove the key and close the door.

**Note:** The outer door can only be locked with the unlock key while the outer door is open. Attempting to turn the unlock key while the outer door is closed may damage the electric lock system.

## Connecting a CO<sub>2</sub> gas cylinder


Before connecting a CO<sub>2</sub> gas cylinder, read the precautions for safe operation on pages 5–9.

1. Get a CO<sub>2</sub> gas cylinder ready and install an optional gas regulator MCO-010R.

**Notes:**

- Use a liquefied-CO<sub>2</sub> gas cylinder (at least 99.5 % pure). The siphon (dip tube) type cannot be used.
- When MCO-010R is not available, install a gas regulator rated at 25 MPa(G) (250 kgf/cm<sup>2</sup>(G), 3600 psi(G)) for the primary side, and 0.25 MPa(G) (2.5 kgf/cm<sup>2</sup>(G), 36 psi(G)) for the secondary side.

2. Using the gas tube provided, connect the connecting port A for the CO<sub>2</sub> gas pipe and the gas regulator of the CO<sub>2</sub> gas cylinder.

**Notes:**

- By using two tube bands provided, connect the gas tube tightly to prevent it from coming off.
- Make sure that the gas tube is not folded.
- If the CO<sub>2</sub> gas is supplied to multiple CO<sub>2</sub> incubators from a single gas cylinder, a CO<sub>2</sub> solid will be formed in the gas regulator. The gas regulator safety valve will operate, and there may be an explosive sound.

3. After connecting the gas tube, make sure that no gas is leaking (ex. by using a gas leak detector spray).

4. Set the CO<sub>2</sub> gas on the secondary side to 0.03 MPa(G)–0.1 MPa(G) (0.3 kgf/cm<sup>2</sup>(G)–1 kgf/cm<sup>2</sup>(G), 4.4 psi(G)–14.5 psi(G)) for gas injection. Recommended pressure: 0.03 MPa(G) (0.3 kgf/cm<sup>2</sup>(G), 4.4 psi(G)).

**Note:** As the pressure increases, the CO<sub>2</sub> gas density variation will increase. Excessive pressure may cause gas supply lines inside the incubator to come loose, which may result in gas poisoning or oxygen deprivation due to leaking gas. If gas lines come loose, the incubator must be repaired.

5. When there is no CO<sub>2</sub> gas left and the CO<sub>2</sub> gas empty alarm is activated, replace the empty gas cylinder.

**Note:** When the optional gas auto-changer MCO-21GC is installed, it switches the empty CO<sub>2</sub> gas supply line to the other automatically. Refer to pages 60–61.

**Note:** The gas lines connected to the incubator will degrade over time. If any deterioration or abnormalities are found during inspection, replace the lines immediately.

# BEFORE COMMENCING OPERATION

## Initial cleaning method

Before using the incubator for the first time, clean dirt (tape residue, smear, etc.) from the chamber and the inner attachments thoroughly.

Keeping the chamber clean is essential for efficient performance. Clean the chamber and all the inner attachments at least once a month. Use the following steps to clean the incubator properly. When these steps are not sufficient to remove dirt, contact our sales representative or agent.



1. Remove the inner attachments, referring to "Removing inner attachments" on pages 23–24.
2. Clean the inner attachments, the chamber inside walls, and the inner door gaskets with a cloth or sponge soaked in neutral detergent, diluted by 5 % or less (Fig. 1). (Undiluted detergent can damage the plastic components. For the dilution, refer to the instructions on the detergent.)

### CAUTION

Do not use detergents or antiseptic solutions with acid, alkali, or chlorine. Doing so may cause discoloration, corrosion, or rusting.

Be careful to keep the diluted detergent or water out of the temperature sensor, the CO<sub>2</sub> gas injection port, the inner sample air access port, the fan motor shaft bearing, and the inner sample air outlet (marked on Fig. 2 by ). Also, do not wash the temperature sensor and the UV lamp using detergent. Doing so may cause failure.

3. Soak a piece of gauze or unwoven cloth in distilled water, wring it out and wipe off the residual detergent thoroughly.
4. Wash the two silicon caps for the access port and the fan using the detergent and rinse them with distilled water, and then autoclave them for sterilization (121 °C, 20 minutes).
5. Wipe the inside walls and the inner attachments thoroughly with an unwoven cloth soaked in alcohol for disinfection. Be careful not to leave any alcohol residue.
6. Reinstall the inner attachments correctly and securely, referring to "Installing inner attachments" on page 25.


## Removing inner attachments


Always replace the inner attachments after cleaning to maintain efficient performance.

Be careful not to damage the UV lamp in the duct (MCO-170AICUVH/170AICUV or when an optional UV system set MCO-170UVS is installed).


1. Turn OFF the power to the incubator.
2. Open the outer and inner doors and pull out all the trays (Fig. 1).


3. Pull out the humidifying pan (Fig. 2).

4. Pull out the humidifying pan cover (Fig. 3).


5. Pull out the fan cover (Fig. 4).


# BEFORE COMMENCING OPERATION



6. Loosen the screws at the rear and lift out the duct (Fig. 5).

7. Pull out the fan (Fig. 6).

8. Remove the silicon caps for the access port from the interior (Fig. 7) and exterior (Fig. 8).

## Installing inner attachments

To reinstall all the attachments, perform the procedure in reverse order from step 8 on page 24.

**Note:** When installing the fan, attach it to the motor shaft securely. Lightly turn and pull the fan manually to make sure that it does not touch the rear panel and is secure (Fig. 1).

### ⚠ CAUTION

If the fan is not pushed fully in, the intended velocity cannot be achieved, which may cause culture failure or insufficient decontamination.

**Note:** To install the duct, check that four screws are securely inserted in the four holes of the duct (Fig. 2).

### ⚠ CAUTION

If the duct is fixed insecurely, the intended velocity cannot be achieved, which may cause culture failure or insufficient decontamination.

**Note:** When installing the fan cover, align the long hole of duct with the projection of the fan cover and push until it clicks into place (Fig. 3). The same applies for the humidifying pan cover.

Make sure that there is no space at the bottom of the fan cover after installing it. If the fan cover is incorrectly installed, it may adversely affect the chamber temperature distribution.

### ⚠ CAUTION

If the fan cover is fixed insecurely, the intended velocity cannot be achieved, which may cause culture failure.

**Note:** Position the tray with the front edge bent down (Fig. 4).

- ① Position the center hole of the fan with the projection of the motor shaft. And insert it deeply.

- ② Lightly turn the fan manually to make sure that it does not touch the rear panel.

- ③ Lightly pull the fan manually to make sure that it is installed.

Fig. 1



Fig. 2


- ① Position and insert.

- ② Confirm the direction of fixing position as shown.

- ③ Make sure that there is no space at the lower of fan cover.

Fig. 3

Bent down

Fig. 4

# BEFORE COMMENCING OPERATION

## Filling the humidifying pan

Use the following procedure to fill the humidifying pan with water or to replace water in it.

1. Pull out the humidifying pan (Fig. 1).
2. Dispose of the remaining water and clean the humidifying pan with a diluted detergent. Then rinse it thoroughly with distilled water and wipe it with alcohol for disinfection.
3. Wipe all moisture from the bottom of the chamber.
4. Return the humidifying pan to the chamber and fill with sterile distilled water (approx. 1.5 L, preheated to 37 °C) (Fig. 2).


Fig. 1


Fig. 2

### Notes:

- Operation with no water for humidifying may increase the chamber temperature than the set temperature temporarily.
- Preheat the water to 37 °C. Adding cool water will lower the temperature and humidity in the chamber.
- Install the humidifying pan in a longitudinal direction with its shorter side placed at the back.
- Refill the humidifying pan with water promptly when the volume of water decreases.
- Mixing any reagent into the water for humidifying may have an adverse effect on the cultivation. Do not add any reagent to the water when using the UV lamp, as the UV light may cause the mixture to deteriorate.

5. Place the humidifying pan with the inner side flush against the back, and close the inner and outer doors.

**Note:** The humidity control bar in the duct maintains low temperature and inner moisture is recondensed. Slide the humidifying pan down right under the humidity control bar, otherwise the recondensed water droplets will fall directly to the bottom of the chamber forming a pool. When the pooled water evaporates, it may leave a white mark on the bottom. This is not a malfunction. Wipe it off with a piece of gauze or unwoven cloth soaked in alcohol. If the mark cannot be removed, scrub it off using a cleansing cream.

# FOR OPTIMAL CULTIVATION

## Precautions for cultures

- **Leave space between culture containers.**

Always leave space for ventilation between culture containers (Petri dishes, flasks, etc.). Inadequate spacing may result in uneven temperature distribution and CO<sub>2</sub> gas density.

- **Do not place harmful materials in the chamber.**

Never place samples that release acidic, alkaline, or corrosive gas into the chamber. Doing so may cause damage resulting from discoloration or corrosion.

- **Close the inner door.**

Always close the inner door before closing the outer door. Failure to close the inner door will adversely affect performance even if the outer door is closed.

- **Open and close the doors gently.**

Closing the doors forcefully may cause spillage of the culture medium, incomplete closing, or damage to the gasket. Before opening the inner door, check through the glass to confirm that the UV lamp is OFF (if the MCO-170AICUVH/170AICUV or the optional MCO-170UVS is installed).

- **Be careful when closing the outer door.**

Use the handle when closing the outer door. Not doing so may cause injury by getting fingers caught in the door. Do not lean on the outer door. Doing so may result in injury from the outer door coming loose or the incubator falling over, or may cause current leakage or electric shock.

- **Be careful of the inside of the outer door.**

The inside of the outer door may become hot.

- **Avoid using excessive force on the inner door.**

Do not put your hand on the glass, poke it with sharp objects, or apply strong force. Doing so may result in injury from breaking the glass.

- **Check the cause of any alarm buzzer.**

If an alarm buzzer sounds while the incubator is in use, check the cause of the alarm immediately. For details on what may cause an alarm buzzer to sound, refer to pages 65–67.

- **Vibration of a shaker.**

When incubators are stacked, operation of a shaker in the chamber of the CO<sub>2</sub> incubator may adversely affect the other incubator.

# FOR OPTIMAL CULTIVATION

## Preventing contamination

To prevent contamination of the chamber, select a suitable installation site.


- **Avoid locations with high temperatures or humidity** where the air may contain more microorganisms.
- **Do not place the incubator near doors, air conditioners, fans, etc.**, where draughts and passage of bodies may increase the risk of microorganisms entering the chamber.
- **If possible, use a cleanroom.**
- **Use clean containers.** Dirty containers are the greatest cause of contamination for cultures. Be careful not to get containers or trays dirty when taking them in and out.
- **Keep the chamber clean.** Wipe off any fingerprints. If water spills from the humidifying pan, or if the doors are left open for a long time, condensation may form on the inside of the doors. Wipe off condensation with a piece of dry sterile gauze. In particular, clean and disinfect the chamber if the culture medium is spilt.
- **Use sterile distilled water in the humidifying pan.** Do not use ultrapure water, which may cause red rust-like particles in the humidifying pan. Clean the humidifying pan once a month.
- **Keep the incubator out of direct airflows from air conditioners or fans.** Cool airflow from an air conditioner may cause condensation and lead to possible contamination.

# CORRECT OPERATION

Use the following procedure to start trial operation or actual operation of the incubator.

1. Install the incubator correctly, referring to pages 17–21.
2. Remove the packing materials from the chamber and inner attachments. Clean and disinfect the chamber and all the inner attachments, referring to page 22.
3. Add approximately 1.5 L of sterile distilled water to the humidifying pan. (Refer to page 26.)
4. Connect the removable power supply cord to the port on the lower right side.

**Note:** Two removable power supply cords are provided.


5. Connect the removable power supply cord to the outlet.
  6. Turn ON the power switch on the lower right side of the incubator.
  7. (For model MCO-170AICUVH/170AICUV or when an optional UV system set MCO-170UVS is installed) Set the frequency of power supply on the LCD touch panel (Refer to page 49).
- The removable power supply cord provided is only for this product. Never use it for any other products.

# BASIC OPERATION ON LCD TOUCH PANEL

## Operation from Menu key

Pressing the Menu key at the bottom right corner of the Top screen brings up the **Menu** screen.


From the Menu screen, the following options are available:

- Press “Set” to bring up the **Stand-by Setting** Screen. Use this to set the temperature, CO<sub>2</sub> density, and high-limit temperature alarm (page 34).
- Press “Log” to bring up the **Log** screen. Use this to display log graphs, export logs, set the log interval, and set the unique ID (pages 38–45).
- Press “Lock” to bring up the **Lock** screen. Use this to set the key lock (page 35), user-ID, password, and optional Auto-lock (pages 56-59).
- Press “Tools#1” to bring up the **Tools#1** screen. Use this to set alarms (pages 36–37), set and operate the optional UV lamp (pages 48–50), and calibrate CO<sub>2</sub> density (pages 63–64).
- Press “Tools#2” to bring up the **Tools#2** screen. Use this to set date, time, brightness, and sleep mode (pages 46–47).

### Notes:

- The **Service** key is available only to qualified engineers.
- Other functions, such as gas supply line changeover and H<sub>2</sub>O<sub>2</sub> decontamination, are performed directly from the Top screen.

To navigate away from any of the screens accessed from the Menu screen, press “Top” to return to the Top screen or “Back” to return to the previous screen.


**Note:**

On the Tools #1 screen, mistakenly pressing the Temp./Gas Calibration key causes the Temp/Gas Calibration screen to be displayed.

When this screen is displayed, use “Back” or “Top” to exit.

Changing the settings on this screen may cause inaccurate temperature or CO<sub>2</sub> density to be displayed.


# BASIC OPERATION ON LCD TOUCH PANEL


## Entering numeric values into input window

On each screen in the LCD touch panel, numerical values may need to be entered on a numerical input box.

1. Pressing a numerical input box causes the input window to be displayed.

2. Press a numerical key or use the Up/Down key (▲/▼) to enter a numerical value. Press "OK".

**Note:** In certain cases, the Up/Down key is not displayed.


**Note:** While the input window is open, it is not possible to operate the Top and Back keys.

## Operating slide keys


On some screens in the LCD touch panel, functions can be turned on and off using a slide key.

To turn a function off, hold the slide key and slide to the left to change the display to OFF. Press “Apply”.

The screenshot shows a dark-themed settings menu titled "Key Lock". At the top right are two buttons: "↑ Top" and "◀ Back". Below the title, there are three input fields: "Key Lock" with a slider set to "OFF", "Password" with an empty text box, and "Confirm Password" with an empty text box. An "Apply" button is located at the bottom right.

To turn a function on, hold the slide key and slide to the right to change the display to ON. Press “Apply”.

The screenshot shows the same "Key Lock" settings menu. The "Key Lock" slider is now set to "ON". The "Password" and "Confirm Password" fields are filled with asterisks (\*\*\*\*). The "Apply" button remains at the bottom right.

# BASIC PARAMETERS

## Setting temperature, CO<sub>2</sub> density, and high-limit temperature alarm

Use the following procedure to set the chamber temperature, the CO<sub>2</sub> density, and the temperature of the high-limit temperature alarm for normal operation. The incubator uses these settings automatically after the power is switched on.

Go to the Stand-by Setting screen (**Menu** → **Set**) and enter each parameter. Press “Apply” to save the entered values.

The screenshot shows the 'Stand-by Setting' screen with the following parameters:

- Temperature : 37.0 °C
- CO<sub>2</sub> Density : 5.0 % (ON)
- High Limit : 53.0 °C

Buttons: ↑ Top, ◀ Back, Apply

### ● Parameter settings

- Temperature: Set value of chamber temperature. Settable range: 0.0 °C–50.0 °C; factory setting: 37.0 °C.
- CO<sub>2</sub> Density: Set value of chamber CO<sub>2</sub> density. Settable range: 0.0 %–20.0 %; factory setting: 0.0 % (OFF). If not necessary to set, slide to OFF, or enter 0.0.
- High Limit: The high-limit temperature alarm is different from the automatic set temperature alarm (page 36), and is independent of it. It is activated if the chamber temperature exceeds the temperature set. Settable range: 20.0 °C–53.0 °C; factory setting: 53.0 °C.

### Notes:

- When changing the set temperature from less than 45.0 °C to 45.0 °C or higher, the incubator readjusts the CO<sub>2</sub> sensor. During readjusting, “Status: Gas sensor initializing” is displayed in the Message display field, and “--.” is displayed in the Present CO<sub>2</sub> density display field. After 15 minutes in the shortest, the incubator returns to the normal operation.
- When operating the incubator for the first time, or after not using it for an extended period, set the chamber temperature to the desired temperature and set the CO<sub>2</sub> density to 0 %. Operate the incubator for at least 4 hours until the chamber temperature and the CO<sub>2</sub> sensor are stable. Then change the setting to the desired CO<sub>2</sub> density.
- Set the high-limit temperature alarm parameter when the chamber temperature is stable at the set value.
- Set the high-limit temperature alarm to at least 1 °C higher than the temperature set for the chamber.

## Setting key lock

When the key lock is set to ON, users are required to enter a password except in the following cases:

- When the CO<sub>2</sub> gas supply line select key on the Top screen is pressed.
- When the Buzzer key on the Top screen is pressed.
- When the optional auto-lock function is ON, to unlock the outer door (see page 58).

Go to the Key Lock screen (**Menu** → **Lock** → **KeyLock**). Set/reset and confirm the password for the key lock. Press “Apply” to turn the key lock ON and to save the password.

Key Lock screen showing the Key Lock toggle set to ON, Password field with four asterisks, and Confirm Password field with four asterisks. The Apply button is visible at the bottom right.

### Notes:

- The maximum length of the password is 6 digits.
- Manage the password for the key lock properly.

## Turning off key lock

Go to the Key Lock screen (**Menu** → **Lock** → **KeyLock**). Slide the Key Lock to OFF and press “Apply” to turn the key lock OFF.

Key Lock screen showing the Key Lock toggle set to OFF, Password field empty, and Confirm Password field empty. The Apply button is visible at the bottom right.

# ALARM PARAMETERS

The incubator is equipped with several alarm functions.

1. Go to the Tools #1 screen (**Menu** → **Tools#1**). Press “Alarm Setting #1” to display the Alarm Setting #1 screen. Set automatic temperature alarm, automatic CO<sub>2</sub> density alarm, and each alarm delay from this screen (procedure 2).

Press “Alarm Setting #2” to display the Alarm Setting #2 screen. It is possible to set door alarm delay, ring back, and remote alarm from this screen (procedure 3).

2. Enter each parameter on the Alarm Setting #1 screen. Press “Apply” to save the entered value. The display returns to the Tools #1 screen.

- Parameter settings

- Temp. Alarm:

Activated when the chamber temperature varies from its set value entered here.

Settable range:  $\pm 1.0$  °C– $\pm 5.0$  °C; factory setting:  $\pm 1.0$  °C.

- CO<sub>2</sub> Alarm:

Activated when the chamber CO<sub>2</sub> density varies from its set value entered here.


Settable range:  $\pm 0.5$  %– $\pm 5.0$  %; factory setting:  $\pm 1.0$  %.

- Alarm Delay:


When the conditions for the automatic temperature alarm or the automatic CO<sub>2</sub> density alarm are reached, the alarm buzzer will sound after the time set here has elapsed.

Settable range: 0 minutes–15 minutes; factory setting: 15 minutes.

**Note:** When the conditions return to the normal range within the alarm delay time, the buzzer does not sound.

3. Enter the required settings for each alarm on the Alarm Setting #2 screen. Press “Apply” to save the settings.

- Settings

- Door Delay:

When the conditions for the door alarm are reached, the alarm buzzer will sound after the time set here has elapsed. Settable range: 1 minute–30 minutes; factory setting: 2 minutes.

**Note:** When the conditions return to normal within the door alarm delay time, the buzzer will not sound.

- Ring Back:

If the alarm buzzer is stopped by pressing the “Buzzer” key, the buzzer will sound again when the conditions that activated the alarm continue after the time set here has elapsed.

Settable range: 1 minute–99 minutes; factory setting: 30 minutes.

**Note:** At Err01 (CO<sub>2</sub> gas cylinder empty), Err11 and Err12 (CO<sub>2</sub> sensor error), Err18 (UV lamp failure), and Door alarm, the alarm is deactivated by pressing the “Buzzer” key and so the buzzer will not sound again (Refer to pages 65–66).

- Remote Alarm:


The remote alarm continues even if the buzzer is stopped by pressing the “Buzzer” key.

Factory setting: ON.

When the alarm sounds, the buzzer can be silenced by pressing “Buzzer”. Refer to Tables 5–7 on page 66 for the way the incubator behaves at the time of pressing “Buzzer” and the re-activation of the alarm under each condition.

Resolve the cause of the alarm, referring to pages 65–67. The alarm itself is not deactivated by pressing “Buzzer” (except for those cases mentioned in the previous note).

**Note:** The buzzer for the high-limit temperature alarm cannot be silenced.

# OPERATION/ALARM LOG

## Setting log interval

The incubator can save operation log data (chamber temperature, CO<sub>2</sub> density, and open/closed state of outer door). Use the following procedure to set the log interval (interval of acquiring the operation log).

1. Go to the Log screen (**Menu** → **Log**). Press “Setting” to display the Setting screen.

2. Enter the log interval. Press “Apply” to save the entered value.

Settable range: 2 minutes–30 minutes; factory setting: 6 minutes.

**Note:** It is possible to register alphanumeric characters (maximum 8 digits) as the Unique ID. Refer to page 42.

**Note:** Relation between log interval and the estimated amount of data that can be saved:

Log interval=2 minutes: Approx. 46 days

Log interval=6 minutes: Approx. 135 days

Log interval=30 minutes: Approx. 664 days

When the amount of data saved exceeds the maximum, the oldest data is deleted and overwritten.


## Displaying operation log

Operation log data can be displayed graphically on the LCD touch panel.

1. Go to the Log screen (**Menu** → **Log**). Press “Chart” to display the Chart screen.


2. Enter the date (year/month/day) of the operation log to be displayed.


3. Press the variable to be displayed and press “Show”. The graph of the operation log is displayed.

- Actual Temp.: (procedure 4)
- Actual CO<sub>2</sub> Level: (procedure 5)
- Door Opening: (procedure 6)


# OPERATION/ALARM LOG


4. Actual Temp. log graph is displayed.

5. Actual CO<sub>2</sub> Level log graph is displayed.

6. Door Opening log graph is displayed.

**Note:** Accurate to about 1 minute per month.


## Exporting operation log

From the log graph screens in the previous page, operation log data can be exported in CSV format to a USB memory device.

1. Insert the USB memory device into the USB port.
2. Press “Export”.

See below for details about abnormal export or name of exported file.

Operation log data can also be exported to a USB memory device from the **Top** Screen.

1. Insert the USB memory device into the USB port.

**Note:** It is not possible to use a USB memory device with a security function that requires entering a password.

2. Go to the Log screen (**Menu** → **Log**). Press “Data Export” to display the Export screen.


3. Select the time period to be exported.

- To export the operation log data for the entire period, press “All” button.

- To export the operation log data for a specified date, press “1 Day” button and enter the date (year/month/day).

**Note:** Accurate to about 1 minute per month.

# OPERATION/ALARM LOG

4. Select the type of log data to be exported.


- To export all types of data, press “All Ch”.
- To export only selected data, press the variable to be exported:
  - Actual Temp.
  - Actual CO<sub>2</sub> Level
  - Door Opening\*

Then press “Selected Ch”.


\* When both of the Auto-lock function and the User-ID mode are ON (refer to pages 56–58), the User-IDs entered to unlock the outer door are also exported.

**Note:** When the specified operation log data does not exist, a warning notice is displayed. Press “OK”, and select a new period, following procedures 3 and 4.

5. When export is completed, press “OK” and remove the USB memory device from the USB port.

**Note:** Even after export is completed, log data saved in the incubator are not deleted.

## Notes:

- A log folder is created in the USB memory device, and the exported file is saved in this folder in CSV format. The name of the exported file is in the format **date\_type of data** (eg 20161001-20170101\_AllCh.csv or 20170101\_Temp.csv).
- On the beginning of the exported file, “MCO-170AIC,” and Unique ID (if registered, refer to page 38) are written (eg MCO-170AIC, RoomA001).

## Displaying alarm log

The incubator can save alarm log data (maximum 256 logs). Alarm log data can be displayed graphically on the LCD touch panel.

**Note:** When number of saved alarm logs exceeds 256, the oldest log is deleted and overwritten.

1. Go to the Log screen (**Menu** → **Log**). Press “Alarm” to display the Alarm screen.

2. Alarm logs from the last 7 days (including current day) are displayed.

**Note:** When the number of alarm logs exceed 6, pressing the top (▲) or the bottom (▼) log, scrolls the log table and allows hidden alarm logs to be seen.

First	Last	Error Code	Warning / Error
2016/12/28 09:38	2016/12/28 09:51		Low CO <sub>2</sub> Density. ▲
2016/12/28 09:25	2016/12/28 09:48	01	CO <sub>2</sub> Gas Empty.
2016/12/28 03:07	2016/12/28 03:55		High Temp.
2016/12/27 22:48	2016/12/27 23:10		High Temp.
2016/12/26 12:19	2016/12/26 12:35		Low CO <sub>2</sub> Density.
2016/12/26 12:02	2016/12/26 12:30	01	CO <sub>2</sub> Gas Empty. ▼

12 / 12

3. Change the number in the “Last X Days” input box to display alarm logs for specified days (including current day).

Settable range: 1 day–45 days.

**Note:** Accurate to about 1 minute per month.

First	Last	Error Code	Warning / Error
2017/01/01 11:31	2017/01/01 11:44		High Temp.

1 / 12


# OPERATION/ALARM LOG

## Exporting alarm log

From the **Alarm** screen, alarm log data can be exported in CSV format to a USB memory device.

1. Insert the USB memory device into the USB port.
2. Press “Export”.

See below for details about abnormal export or name of exported file.

Alarm log data can also be exported to a USB memory device from the **Top** Screen.


1. Insert the USB memory device in the USB port.
2. Go to the Log screen (**Menu** → **Log**). Press “Alarm Export” to display the Alarm Export screen.


3. Select the period for export.

- For entire period, press “All”.
- For specified recent days (including current day), press “Last X Days” and enter number.  
Settable range: 1 day–45 days.


**Note:** Accurate to about 1 minute per month.

4. Press “Export”.


**Note:** A warning notice appears when no alarm log data exists for the specified days. Press “OK” and enter a new period, following procedure 3.

5. When export is completed, press “OK” and remove the USB memory device from the USB port.

**Note:** Even after export is completed, log data saved in the incubator will not be deleted.

**Note:** A log folder is created in the USB memory device, and the exported data file is saved in this folder in CSV format. The name of the exported file is in the format **earliest date-last date\_AlarmLog**. (eg 20170101-20170107\_AlarmLog.csv).

# OTHER PARAMETERS

## Setting date and time


1. Go to the Tools #2 screen (**Menu** → **Tools#2**). Press “Date & Time” to display the Date & Time screen.

2. Enter the present date and time. Press “Apply” to save the entered value.

**Notes:**

- 24-hour clock.
- Accurate to about 1 minute per month. It is recommended that the time is reset periodically.


## Setting brightness and sleep

1. Go to the Tools #2 screen (**Menu** → **Tools#2**). Press “Brightness/Sleep” to display the Brightness/Sleep screen.

2. Select the required settings for brightness and sleep. Press “Apply” to save the settings.

- Settings

- Brightness(Active):

Brightness of LCD touch panel in the normal state. Adjust the slide bar or enter a numerical value into the input box. Settable range: 50–100; factory setting: 80.

- Sleep:

Lowers brightness of the LCD touch panel to save electricity during periods of inactivity.

Sleep function is turned ON or OFF using the slide key. Enter a value to set the time before sleep mode occurs. Settable range: 1 minute–5 minutes; factory setting: 2 minutes.

**Note:** It is not possible to operate any key when the LCD panel is in sleep mode. Touching the panel releases the sleep mode and allows normal operation.

- Brightness(Sleep):

Brightness of LCD touch panel during sleep mode. Adjust the slide bar or enter a numerical value into the input box. Settable range: 0–50; factory setting: 20.

# UV LAMP PARAMETERS

The UV lamp operates in models MCO-170AICUVH/170AICUV or when the optional UV system MCO-170UVS is installed on the model MCO-170AIC.

Operation of the UV lamp disinfects the water in the humidifying pan and the air circulating in the chamber.

## Using the UV lamp

1. Correctly install all of the inner attachments, and place the cultivation samples on the trays.

**Notes:**

- The duct, humidifying pan, and humidifying pan cover prevent UV light from escaping. Always use the humidifying pan and cover, even when not humidifying.
- Never turn the UV lamp ON when the duct and humidifying pan cover are removed.
- Always use the humidifying pan cover even when the UV lamp is not turned ON. Using the incubator without the humidifying pan cover may adversely affect the chamber temperature distribution and humidity recovery.

2. Close the outer door. The UV lamp will light for the preset period (i.e. the time set on “UV Timer” plus the time added by “UV Timer Ext”). Refer to page 49 for instructions on setting the timer.

**Notes:**

- If the outer door is opened while the UV lamp is lit, the lamp will turn OFF. When the door is closed again, the lamp will light for the preset period.
- Repeated opening and closing of the outer door may adversely affect the condensation in the chamber and the chamber temperature distribution because the UV lamp generates heat for a long time. It may also shorten the service life of the UV lamp.
- To check whether the UV lamp is lit, open the outer door and then press the door switch with the inner door closed. Visible blue light can be observed from the front of the humidifying pan cover.
- UV light is harmful to the eyes. Do not turn the UV lamp ON when the inner door is open.

3. If the outer door is not opened for at least 12 consecutive hours, the UV lamp lights for the preset period every 12 hours.

**Note:** Opening the outer door will reset the 12-hour period.

- The recommended replacement time for the UV lamp (i.e., when the UV output ratio drops to 60 %–70 % of its initial value) is when the accumulated ON time reaches 5,000 hours. The message “Warning: UV Bulb Life” then appears in the message display field. Replace the UV lamp promptly at this point. Contact our sales representative or agent for information.


- If the UV lamp burns out, “Err18: UV Lamp Abnormal” is displayed in the message display field. If this occurs, it is not possible to perform H<sub>2</sub>O<sub>2</sub> decontamination. Replace the UV lamp, and replace the glow starter (type FG-7P or FG-7PL) at the same time. Contact our sales representative or agent for information.


## Setting UV lamp ON period

Use the following procedure to change the setting of the UV lamp ON period:

1. Go to the Tools #1 screen (**Menu** → **Tools#1**). Press “UV Setting” to display the UV Setting screen.

2. Enter the required settings. Press “Apply” to save the settings.

- Settings
- UV Timer:

The time for which the UV lamp remains lit after closing the outer door.

Settable range: 0 minutes–30 minutes; factory setting: 10 minutes.

**Notes:**

- It is recommended that the UV Timer is set for 10 minutes. A lower setting may result in insufficient disinfection.
- When the UV timer is set at 0, the UV lamp does not light.

- UV Life:

The total time for which the UV lamp has been lit is displayed as a percentage of the 5,000 hours before replacement is recommended (It is impossible to alter this setting).

- UV Timer Ext.:

As the UV lamp is used, its UV ray output progressively declines. To compensate for this, the lighting time is automatically extended as the lifetime of the UV lamp increases (The setting of UV Timer is unchanged). Extension rate: 0 %–40 % (impossible to alter); factory setting: 0 %.

For example, if UV Timer is set to 10 minutes and UV Timer Ext. is 40 %, then the UV lamp lights for 14 minutes.

- Frequency:

Frequency of the power supply. Press button for 50 Hz or 60 Hz. Factory setting: 50 Hz.

# UV LAMP PARAMETERS

## Lighting the UV lamp for 24 hours


If the chamber has been contaminated by dirt or by spilling the medium, use the following procedure to decontaminate it by lighting the UV lamp for 24 hours.

1. Remove all attachments from the chamber, including the trays, the fan cover, the duct, the fan, the humidifying pan, and the humidifying pan cover. Disinfect all the attachments in an autoclave or with alcohol.
2. Clean and wipe down the inside the chamber with alcohol.
3. Set the CO<sub>2</sub> density to 0 %. Refer to page 34.

4. Go to the Tools #1 screen (**Menu** → **Tools#1**). Press “UV Setting” to display the UV Setting screen.

5. Turn the UV 24h Mode to ON and press “Apply”.

6. The UV lamp will now light continuously for 24 hours. “UV : ON” is displayed on the UV lamp condition display.

### Notes:

- The UV 24-hour mode may activate the automatic set temperature alarm because of rising chamber temperature.
- If the outer door is opened when the UV lamp is lit, the UV lamp is turned OFF and UV 24-hour mode is cancelled. Repeat procedures 4 to 6 to restart the UV 24-hour mode.

7. After 24 hours, the UV lamp turns OFF automatically. Reinstall all the attachments removed in procedure 1.

# H<sub>2</sub>O<sub>2</sub> DECONTAMINATION

H<sub>2</sub>O<sub>2</sub> decontamination can be performed when the chamber is contaminated or when cleaning the chamber prior to starting a culture.

For H<sub>2</sub>O<sub>2</sub> decontamination to be performed, one of the following conditions must apply:

- An H<sub>2</sub>O<sub>2</sub> generator MCO-HP is installed in the MCO-170AICUVH.
- An H<sub>2</sub>O<sub>2</sub> generator MCO-HP, H<sub>2</sub>O<sub>2</sub> decon board MCO-170HB, and electric lock MCO-170EL are all installed on the model MCO-170AICUV.
- UV system set MCO-170UVS, H<sub>2</sub>O<sub>2</sub> generator MCO-HP, H<sub>2</sub>O<sub>2</sub> decon board MCO-170HB, and electric lock MCO-170EL are all installed on the model MCO-170AIC.

When optional small doors (MCO-170ID) are installed, H<sub>2</sub>O<sub>2</sub> decontamination can be performed with an MCO-170ID.

Read the precautions for safe operation on pages 5–9 before performing decontamination.

## Notes:


- Never start H<sub>2</sub>O<sub>2</sub> decontamination when “Warning: UV Bulb Life” is displayed in the message display field. The UV resolution will be insufficient.
- If the power supply is cut off during H<sub>2</sub>O<sub>2</sub> decontamination, the electric lock will remain ON. After the power supply is restored, the H<sub>2</sub>O<sub>2</sub> gas resolution process will restart and finish automatically. In this case, repeat the decontamination process because the decontamination will not be complete.

## H<sub>2</sub>O<sub>2</sub> decontamination

1. Take all the trays, the fan cover, the duct, the humidifying pan cover, and the humidifying pan out of the chamber.
2. Dispose of the water in the humidifying pan, and wipe the attachments removed from the chamber with a piece of gauze soaked with water or alcohol for disinfection.
3. Wipe the inside walls of the chamber with a piece of gauze soaked with water or alcohol.
4. Attach the duct and the fan cover.
5. Pour one bottle of the H<sub>2</sub>O<sub>2</sub> reagent MCO-H2O2 into the H<sub>2</sub>O<sub>2</sub> generator MCO-HP (Fig. 1).

6. Fit the two pins on the H<sub>2</sub>O<sub>2</sub> generator into the two holes on the lower left side of the duct (Fig. 1).

**Note:** Make sure that the H<sub>2</sub>O<sub>2</sub> generator is securely attached. Improper attachment may result in insufficient decontamination.

# H<sub>2</sub>O<sub>2</sub> DECONTAMINATION

7. Remove the connector cap from the connector on the bottom right of the far side of the chamber. Connect the H<sub>2</sub>O<sub>2</sub> generator and connector by the cable provided (Fig. 2).

**Note:** Be sure to keep the connector cap.


Fig. 2

8. Insert four trays in the 2nd, 3rd, 4th, and 6th tray catches from the top of the chamber side.

**Note:** The trays included as accessories are designed to be appropriate for decontamination. If half trays (MCO-25ST, optional) or trays for previous models are used, decontamination may not be sufficiently effective.

9. Insert the humidifying pan cover and humidifying pan removed in procedure 1 (Fig. 3).



Fig. 3


10. Make sure that the duct, the fan, and the H<sub>2</sub>O<sub>2</sub> generator are securely attached (refer to page 25). Then close the inner and outer doors.

**Note:** When optional small doors (MCO-170ID) are installed, make sure that the 4 gastight split doors are closed.


11. Press “H<sub>2</sub>O<sub>2</sub>” for 3 seconds to open the Setting Position window. Enter the password if prompted.

12. Press “Next” to display the H<sub>2</sub>O<sub>2</sub> Decontamination Step1 screen. If there is something to worry about setting position, check the inside again.


13. The system check starts automatically. If an alarm or safety function is activated, refer to Table 8 on page 67.

14. If the system is normal, the H<sub>2</sub>O<sub>2</sub> Decontamination Step2 screen appears. Press “OK” to start H<sub>2</sub>O<sub>2</sub> decontamination. H<sub>2</sub>O<sub>2</sub> decontamination is performed automatically from Step3 to Step8 (procedure 15).

**Notes:**


- The outer door is electrically locked for safety until H<sub>2</sub>O<sub>2</sub> decontamination is completed.
- A buzzer sounds when H<sub>2</sub>O<sub>2</sub> decontamination is completed (after about 100 minutes).

# H<sub>2</sub>O<sub>2</sub> DECONTAMINATION


- Step3 is the process of warming the chamber temperature to 45 °C.

**Note:** By pressing “Abort”, decontamination will be stopped and the procedure jumps to Step8 (end of decontamination).


- Step4 is the process of decontaminating the chamber by generating H<sub>2</sub>O<sub>2</sub> vapour from the H<sub>2</sub>O<sub>2</sub> generator.

**Note:** By pressing “Abort”, decontamination will be stopped and the procedure jumps to Step7 (resolution of H<sub>2</sub>O<sub>2</sub> vapour).


- Step7 is the process of resolving H<sub>2</sub>O<sub>2</sub> vapour using the UV lamp.

**Note:** Step5 and Step6 do not exist.

15. When H<sub>2</sub>O<sub>2</sub> decontamination is completed, the H<sub>2</sub>O<sub>2</sub> Decontamination Step8 screen is displayed. Open the outer and inner doors. Disconnect the connector, cable, and H<sub>2</sub>O<sub>2</sub> generator, and remove the H<sub>2</sub>O<sub>2</sub> generator and cable from the chamber.

**Note:** Put on protective glasses and rubber gloves to do this.

16. Press "OK" to return to the Top screen.

**Note:** The OK key does not work when the H<sub>2</sub>O<sub>2</sub> generator is attached to the connector by the cable.

17. Dilute the remaining H<sub>2</sub>O<sub>2</sub> reagent in the H<sub>2</sub>O<sub>2</sub> generator with a large volume of water and dispose of it. Rinse and wash the H<sub>2</sub>O<sub>2</sub> generator with distilled water. Keep the H<sub>2</sub>O<sub>2</sub> generator in a clean environment outside the chamber.

**Note:** Do not wash either the inside or outside of the H<sub>2</sub>O<sub>2</sub> generator with alcohol.

18. After H<sub>2</sub>O<sub>2</sub> decontamination, surplus H<sub>2</sub>O<sub>2</sub> liquid will remain at the bottom of the chamber and in the bottom part of the H<sub>2</sub>O<sub>2</sub> generator duct. This solution contains H<sub>2</sub>O<sub>2</sub> at a low density. Wear protective glasses and rubber gloves and wipe it up with a non-woven cloth.

19. Ventilate the chamber sufficiently and place all the attachments back inside the chamber.

**Note:** After H<sub>2</sub>O<sub>2</sub> decontamination, replace the connector cap firmly over the connector on the side of the chamber (Fig. 4).


Fig. 4

## Precautions when handling H<sub>2</sub>O<sub>2</sub> reagent

Observe the following precautions when handling optional H<sub>2</sub>O<sub>2</sub> reagent MCO-H2O2:

- Handling precautions

- Wear protective equipment, such as protective glasses and rubber gloves.
- Do not start a fire in the area where the reagent is being handled.
- Do not leave any reagent in the container after it has been used or while it is being used.
- Do not place flammable or combustible materials near the area where the reagent is handled.

- Precautions for storage

- Store in a cool, dark place.
- Always close the container cover securely to prevent impurities from becoming mixed into the reagent.
- Check the container to make sure that there is no damage, corrosion, or cracking.
- Store the container with the inlet facing upwards, and make sure that the container will not tip or be knocked over.

- Precautions for disposal

- Dispose in accordance with local laws and regulations.

# ELECTRIC LOCK (OPTIONAL)

Automatic electric locking is available on model MCO-170AICUVH, or when the optional electric lock MCO-170EL is installed on the models MCO-170AICUV/170AIC.


There are two modes to unlock the outer door.

- Quick mode: Press “Unlock” on the Top screen.
- User-ID mode: Enter the User-ID and password for Auto-Lock, after pressing “Unlock” on the Top screen.

## Setting User-ID

Before turning the User-ID mode to ON, register a User-ID and a password for Auto-Lock (maximum 99 User-IDs and its passwords).


1. Go to the Lock screen (**Menu** → **Lock**). Press “Auto-Lock User” to display the Auto-Lock User screen.


2. Enter desired User-ID and password, and confirm. Press “Add” to register them.

### Notes:

- User-ID: The alphanumeric characters (maximum 8 digits).
- Password: The number (maximum 6 digits).
- A User-ID can be registered without its password.
- Manage the password for Auto-Lock properly.

**Note:** When registering the 100th User-ID, a warning notice is displayed. Press “OK”, and delete a disused User-ID in reference to below.

- How to delete the User-ID already registered

Enter the User-ID to delete, and enter its password in the “Password” input box. Press “Delete”.

**Note:** When deleting all registered User-IDs, the User-ID mode is turned to OFF.


- How to change the User-ID's password already registered

Enter the User-ID whose password to change, and enter its new password in both "Password" and "Confirm Password" input boxes. Press "Add" to re-write the new password.

## Setting auto-lock

1. Go to the Lock screen (**Menu** → **Lock**). Press "Auto-Lock" to display the Auto-Lock screen.

2. Enter the required setting. Then press "Apply".

- Settings
- Auto-Lock:


When Auto-lock is ON, the outer door is locked automatically when the set time has elapsed after the door is closed. Settable range: 1 minute–60 minutes; factory setting: 1 minute.

- User-ID:

Choose the quick mode (OFF) or the User-ID mode (ON). Factory setting: OFF.

### Notes:

- When no User-ID is registered, a warning notice is displayed. Press "OK", and register new User-ID and its password in reference to pages 56–57.
- Even after the User-ID mode is turned to OFF, the User-IDs already registered are not deleted.


# ELECTRIC LOCK (OPTIONAL)

- How to unlock the outer door

- In the quick mode: Press “Unlock” on the Top screen.

**Note:** When the outer door is closed again, it re-locks automatically when the set time has elapsed.


- In the User-ID mode: Press “Unlock” on the Top screen, and enter the User-ID and its password for Auto-Lock.

**Note:** The User-ID entered at this time is saved as the open/close state of outer door log data (refer to pages 41–42).

- When the entered User-ID does not exist, a warning notice is displayed. Press “OK”, and enter the correct User-ID.


## Using unlock key

For model MCO-170AICUVH, or when the optional electric lock MCO-170EL is installed on models MCO-170AICUV/170AIC, the outer door is electrically locked during a power outage. Use the unlock key provided to unlock the outer door when the power is interrupted. To re-lock the outer door, turn the unlock key in the direction of the lock while the outer door is open. Then close the outer door.


**Note:** The outer door can only be locked with the unlock key while the outer door is open. Attempting to turn the unlock key while the outer door is closed may damage the electric lock system.

## Turning off auto-lock

1. Go to the Lock screen (**Menu** → **Lock**). Press “Auto-Lock” to display the Auto-Lock screen.

2. Hold the Auto-Lock slide key and slide it to the left. Then press “Apply”.


# GAS AUTO-CHANGER (OPTIONAL)

When the optional gas auto-changer MCO-21GC is installed, there will be two connecting ports for the CO<sub>2</sub> gas pipe: A and B. When two CO<sub>2</sub> gas cylinders are connected, the auto-changer switches the CO<sub>2</sub> gas supply line when one of the cylinders becomes empty.

## Connecting CO<sub>2</sub> gas cylinders

Refer to page 21 for the procedure for connecting gas cylinders.

1. Take two CO<sub>2</sub> gas cylinders (A and B) and install the optional gas regulator MCO-010R on both.
2. Using one of the gas tubes provided, connect the gas regulator of cylinder A to connecting port A.
3. Using another gas tube, connect the gas regulator of cylinder B to connecting port B.

## Automatic CO<sub>2</sub> gas supply line changeover

When the optional gas auto-changer MCO-21GC is installed, CO<sub>2</sub> gas supply line indicators A and B and the CO<sub>2</sub> gas supply line select key will be displayed on the Top screen. A light will show which CO<sub>2</sub> gas supply line is being used (see page 14).

When the CO<sub>2</sub> density level remains unchanged, even though the CO<sub>2</sub> gas valve in the unit is open, the unit regards the currently connected CO<sub>2</sub> gas cylinder as empty. The CO<sub>2</sub> gas supply line is changed over automatically. The process is set out in Table 2.

1. When there is gas remaining in cylinder A, the unit operates with gas supplied from cylinder A (Situation 1).
2. When cylinder A is empty, the level of CO<sub>2</sub> density in the unit does not increase because gas is not supplied into the unit even though the gas valve in the unit is open (Situation 2).
3. When Situation 2 continues for 2–3 minutes, the gas supply line is changed over automatically. The CO<sub>2</sub> gas empty alarm is activated sounding the buzzer, and gas supply indicator A is displayed in reverse video and blinks (Situation 3).
4. Stop the CO<sub>2</sub> gas empty alarm by pressing “Buzzer”. The light on gas supply indicator A is turned off (Situation 4).

5. Exchange the empty gas cylinder A for a new one immediately (Situation 5).

6. When cylinder B is empty, the supply reverts to cylinder A.

**Table 2. CO<sub>2</sub> gas supply line automatic changeover**

Situation	CO <sub>2</sub> gas			CO <sub>2</sub> gas supply line indicator			Message display field
	Supply line	Cylinder A	Cylinder B		A	B	
<b>1</b> CO <sub>2</sub> gas is supplied from valve A	A	Remaining	Remaining		Light on	Light off	
<b>2</b> CO <sub>2</sub> density in the chamber is not increased even though the CO <sub>2</sub> gas valve is open	A	Empty	Remaining		Light on	Light off	
<b>3</b> CO <sub>2</sub> gas supply line is automatically changed to valve B	B	Empty	Remaining		Reverse video and blink	Light on	Err01: CO <sub>2</sub> Gas Empty (and buzzer)
<b>4</b> "Buzzer" is pressed	B	Empty	Remaining		Light off	Light on	
<b>5</b> New cylinder A fitted	B	Remaining	Remaining		Light off	Light on	

**Notes:**


- When the Buzzer key is not pressed in Situation 4 and cylinder B becomes empty before cylinder A is replaced, the supply will switch back to cylinder A. In this case, replace both cylinders and press "Buzzer" immediately.
- The automatic changeover of cylinders is linked to the CO<sub>2</sub> density in the chamber. The automatic changeover of cylinders may also occur in the situations where, e.g. the gas tube gets clogged, gas is leaking, the gas pressure has dropped, the valve on the gas cylinder is insufficiently open etc, even though the cylinder being used is not completely empty.

# GAS AUTO-CHANGER (OPTIONAL)


## Manual CO<sub>2</sub> gas supply line A/B changeover

It is possible to change the CO<sub>2</sub> gas supply line A/B manually at any time. For example:

1. Press the gas supply line select key "B" for a few seconds.

2. CO<sub>2</sub> gas supply line is changed from A to B.

**Note:** Table 3 sets out the process if the gas supply has been changed automatically and supply line is then manually changed to an empty cylinder before the Buzzer key is pressed.

**Table 3. Manual change from cylinder B to empty cylinder A**

	Situation	CO <sub>2</sub> gas			CO <sub>2</sub> gas supply line indicator			Message display field
		Supply line	Cylinder A	Cylinder B		A	B	
1	Supply has been automatically changed to supply line B	B	Empty	Remaining		Reverse video and blink	Light on	Err01: CO <sub>2</sub> Gas Empty (and buzzer)
2	Buzzer key not pressed but gas supply line select key is pressed	A	Empty	Remaining		Blink	Light off	Err01: CO <sub>2</sub> Gas Empty (and buzzer)

# STD GAS AUTO-CALIBRATION KIT (OPTIONAL)

When the STD gas auto-calibration kit MCO-SG is installed, it is possible to calibrate CO<sub>2</sub> density manually.

1. Connect a standard gas cylinder to the connecting port on the lower right side of the incubator. Use a cylinder with the same CO<sub>2</sub> density as the density set on the incubator.

**Note:** There is no problem in keeping the standard gas cylinder connected after CO<sub>2</sub> gas density calibration is completed.



2. Go to the Tools #1 screen (**Menu** → **Tools#1**). Press “STD Gas Setting” to display the STD Gas Setting screen.

3. Enter the CO<sub>2</sub> density of the connected standard CO<sub>2</sub> gas cylinder. Press “Apply”. The display returns to the Tools #1 screen.

Settable range: 4.0 %–21.0 %; factory setting: 5.0 %.


**Note:** Don't turn the CO<sub>2</sub> STD Gas Density to OFF.

4. Press “STD Gas Calibration” to display the STD Gas Calibration screen.


# STD GAS AUTO-CALIBRATION KIT (OPTIONAL)


5. The system check starts automatically. If the system is normal, the display goes from Step1 to Step2.


6. Press "OK" to go to Step3.

7. At Step3, CO<sub>2</sub> gas density calibration starts. The display proceeds automatically to Step5.

8. After CO<sub>2</sub> density calibration finishes, the incubator returns to the normal operation. Press "OK".

# ALARMS, SAFETY, AND SELF-DIAGNOSIS

The incubator supports the following alarms, safety functions, and self-diagnostic functions. If an error from Err05 to Err18 is activated, or if Err56 is activated, contact our sales representative or agent.

**Table 4. Alarms, safety, and self-diagnosis for culture operations**

Alarm or safety function	Conditions	Display	Buzzer	Remote alarm	Safety operation
High-limit temperature alarm	The chamber temperature exceeds the high-limit alarm temperature set value.	"Over Heat" is displayed alternately in normal characters and reverse video in the Over heat display.	Continuous tone	ON	Heater OFF
Automatic set temperature alarm	The chamber temperature is out of the automatic set temperature alarm setting range ( $\pm 1.0$ °C to $\pm 5.0$ °C).	"Warning :High Temp" or "Warning :Low Temp" is displayed in the message display field.	Intermittent tone after set alarm delay time (0 min to 15 min) has elapsed	ON after set alarm delay time (0 min to 15 min) has elapsed	-----
Automatic set CO <sub>2</sub> density alarm	The chamber CO <sub>2</sub> density is out of the automatic set CO <sub>2</sub> density alarm setting range ( $\pm 0.5$ % to $\pm 5.0$ %).	"Warning :High CO <sub>2</sub> Density" or "Warning :Low CO <sub>2</sub> Density" is displayed in the message display field.	Intermittent tone after set alarm delay time (0 min to 15 min) has elapsed	ON after set alarm delay time (0 min to 15 min) has elapsed	-----
Auto-return	On screens other than the Top screen, there is no key operation for approximately 90 seconds. After sleep function is turned ON, there is no alarm/error and key operation for approximately 90 seconds.	(Return to the "Top screen")	-----	-----	-----
Door alarm	The outer door is open.	"Door : Open" is displayed alternately in normal characters and reverse video in the outer door (opening/closing) display.	Intermittent tone after set alarm delay time (1 min to 30 min) has elapsed	-----	The CO <sub>2</sub> valve is closed. The heater turns OFF after 1 min.
Door lock error	(MCO-170AICUVH or when optional MCO-170EL is installed) Outer door is opened when it is auto-locked by electric lock.	"Err20: Door Lock Failure" is displayed in the message display field.	Intermittent tone	ON	UV lamp OFF
CO <sub>2</sub> gas cylinder empty	The CO <sub>2</sub> density does not increase when the CO <sub>2</sub> valve is opened.	"Err01: CO <sub>2</sub> Gas Empty" is displayed in the message display field.	Intermittent tone	ON	-----
Chamber temperature sensor error	The chamber temperature sensor is disconnected.	"Err05: Temp Sensor Open" is displayed in the message display field.	Intermittent tone	ON	Heater OFF
	The chamber temperature sensor is short-circuited.	"Err06: Temp Sensor Short" is displayed in the message display field.	Intermittent tone	ON	Heater OFF
Sensor box temperature sensor error	The sensor box temperature sensor is disconnected.	"Err07: CO <sub>2</sub> Box Temp Sensor Open" is displayed in the message display field.	Intermittent tone	ON	CO <sub>2</sub> valve is closed.
	The sensor box temperature sensor is short-circuited.	"Err08: CO <sub>2</sub> Box Temp Sensor Short" is displayed in the message display field.	Intermittent tone	ON	CO <sub>2</sub> valve is closed.
Ambient temperature sensor error	The ambient temperature sensor is disconnected.	"Err09: AT Sensor Open" is displayed in the message display field.	Intermittent tone	ON	-----
	The ambient temperature sensor is short-circuited.	"Err10: AT Sensor Short" is displayed in the message display field.	Intermittent tone	ON	-----
CO <sub>2</sub> sensor error	The Vref or Vgas output voltage for the CO <sub>2</sub> sensor is abnormal.	"Err11: CO <sub>2</sub> Sensor Vref Abnormal" or "Err12: CO <sub>2</sub> Sensor Vgas Abnormal" is displayed in the message display field.	Intermittent tone	ON	CO <sub>2</sub> valve is closed.
Main heater error	Main heater burnout occurs or the main heater SSR is short-circuited.	"Err13: Main Heater Abnormal" is displayed in the message display field.	Intermittent tone	ON	-----
Bottom heater error	Bottom heater burnout occurs or the bottom heater SSR is short-circuited.	"Err14: Humidity Heater Abnormal" is displayed in the message display field.	Intermittent tone	ON	-----
Door heater error	Door heater burnout occurs or the door heater SSR is short-circuited.	"Err15: Door Heater Abnormal" is displayed in the message display field.	Intermittent tone	ON	-----
Sensor box heater error*1	a) High-limit temperature alarm is activated. b) Sensor box heater burnout occurs or the sensor box heater SSR is short-circuited.	"Err16: CO <sub>2</sub> Box Heater Abnormal" is displayed in the message display field.	Intermittent tone	ON	-----
Heater SSR burnout*1	a) High-limit temperature alarm is activated. b) Main, bottom, door, or sensor box heater SSR burnout occurs.	"Err17: Heater SSR Open" is displayed in the message display field.	Intermittent tone	ON	-----

\*1: Shortly after the high-limit temperature alarm is activated, Err16 (Sensor box heater error) and Err17 (Heater SSR burnout) are activated.

# ALARMS, SAFETY, AND SELF-DIAGNOSIS

Alarm or safety function	Conditions	Display	Buzzer	Remote alarm	Safety operation
UV lamp failure	(MCO-170AICUVH/MCO-170AICUV, or when optional MCO-170UVS is installed) The UV lamp burns out.	"Err18: UV Lamp Abnormal" is displayed in the message display field.	Intermittent tone	ON	----
New UV lamp replacement	(MCO-170AICUVH/MCO-170AICUV, or when optional MCO-170UVS is installed) The accumulated ON time reaches approx. 5,000 hours.	"Warning: UV Bulb Life" is displayed in the message display field.	----	----	----
Communication error	Communication between LCD touch panel and control substrate is terminated or unstable.	"Err56: Communication Failure" is displayed in the message display field.	----	----	----
Warming-up of gas control	After power switch is turned ON, under warming-up before temperature is stable and gas control is enabled.	"Status: Gas sensor initializing" is displayed in the message display field.	----	----	----

- Tables 5–7 show the behaviour of the alarm (buzzer) and ring back function when pressing the Buzzer key.

**Table 5. In cases other than those covered in Tables 6 and 7**

Remote Alarm setting	Ring Back setting	Buzzer from CO <sub>2</sub> incubator		Remote Alarm	
		When pressing the Buzzer key	When the Ring Back set time passes	When pressing the Buzzer key	When the Ring Back set time passes
ON: Remote alarm setting not connected with Buzzer key (remote alarm is cancelled)	ON	OFF (Alarm is not cancelled)	ON	ON	ON (Under continuation)
	OFF		OFF		
OFF: Remote alarm setting connected with Buzzer key (remote alarm is not cancelled)	ON		ON	OFF (Alarm is not cancelled)	ON
	OFF		OFF		OFF

**Note:** Resolve the cause of the alarm, as the alarm itself is not deactivated by pressing the Buzzer key.

**Table 6. In cases of high-limit temperature alarm or Err38**

Remote Alarm setting	Ring Back setting	Buzzer from CO <sub>2</sub> incubator		Remote Alarm	
		When pressing the Buzzer key	When the Ring Back set time passes	When pressing the Buzzer key	When the Ring Back set time passes
ON: Remote alarm setting not connected with Buzzer key (remote alarm is cancelled)	ON	ON	ON (Under continuation)	ON	ON (Under continuation)
	OFF				
OFF: Remote alarm setting connected with Buzzer key (remote alarm is not cancelled)	ON			(Continue)	(Continue)
	OFF				

**Note:** Close the outer door when Err38 is activated.

**Table 7. In cases of Err01, Err11, Err12, Err18, or door alarm\*2**

Remote Alarm setting	Ring Back setting	Buzzer from CO <sub>2</sub> incubator		Remote Alarm	
		When pressing the Buzzer key	When the Ring Back set time passes	When pressing the Buzzer key	When the Ring Back set time passes
ON: Remote alarm setting not connected with Buzzer key (remote alarm is cancelled)	ON	OFF (Alarm is cancelled)	OFF (Alarm is already cancelled)	OFF (Alarm is cancelled*2)	OFF (Alarm is already cancelled*2)
	OFF				
OFF: Remote alarm setting connected with Buzzer key (remote alarm is not cancelled)	ON				
	OFF				

\*2: The remote alarm does not work for the door alarm.

**Note:** When Err01 is activated, connect the new CO<sub>2</sub> gas cylinder and press the Buzzer key to stop the buzzer. In addition, when the optional MCO-21GC is installed and the gas supply is switched to the reserve gas cylinder, press the Buzzer key and replace the gas cylinder.

**Table 8. Alarms and safety functions for H<sub>2</sub>O<sub>2</sub> decontamination**

Alarm or safety function	Conditions	Display	Buzzer	Remote alarm	Safety operation
System check error at start of H <sub>2</sub> O <sub>2</sub> decontamination	The H <sub>2</sub> O <sub>2</sub> generator is not connected.	"Err31:H2O2 Unit Not Connect" is displayed in the message display field.	-----	-----	H <sub>2</sub> O <sub>2</sub> decontamination is cancelled.
	There is no H <sub>2</sub> O <sub>2</sub> solution or the H <sub>2</sub> O <sub>2</sub> level sensor has failed (or is disconnected).	"Err32: Low H2O2 Level" is displayed in the message display field.	-----	-----	H <sub>2</sub> O <sub>2</sub> decontamination is cancelled.
	The door is not closed.	"Err33: Outer Door Open" is displayed in the message display field.	-----	-----	H <sub>2</sub> O <sub>2</sub> decontamination is cancelled.
Failure during warming-up	Power was interrupted.	After power is restored, the Top screen is displayed.	-----	-----	-----
	The outer door is open.	"Err38:Door Lock Failure" is displayed in the message display field.	Continuous tone (when the outer door is open)	ON (when the outer door is open)	H <sub>2</sub> O <sub>2</sub> decontamination is cancelled.
Failure during H <sub>2</sub> O <sub>2</sub> decontamination	The volume of H <sub>2</sub> O <sub>2</sub> generated is abnormal.	"Err34:H2O2 Volume" is displayed in the message display field.	Intermittent tone with 15-minute delay	ON with 15-minute delay	Moves to UV resolution.
	Power was interrupted during H <sub>2</sub> O <sub>2</sub> decontamination.	After power is restored, "Err35: Power Failure" is displayed in the message display field.	Intermittent tone	ON	<ul style="list-style-type: none"> <li>• During power failure, outer door is electrically locked.</li> <li>• After power is restored, moves to UV resolution.</li> </ul>
	The outer door is open.	"Err38:Door Lock Failure" is displayed in the message display field.	Continuous tone (when the outer door is open)	ON (when the outer door is open)	<ul style="list-style-type: none"> <li>• H<sub>2</sub>O<sub>2</sub> mist generation is cancelled.</li> <li>• Resolution time is extended.</li> </ul>
Failure during H <sub>2</sub> O <sub>2</sub> gas resolution by UV	The UV lamp failed during H <sub>2</sub> O <sub>2</sub> gas resolution by UV.	"Err36: UV Lamp Failure" is displayed in the message display field.	Intermittent tone	ON	Resolution time is extended.
	Power was interrupted during H <sub>2</sub> O <sub>2</sub> gas resolution by UV.	After power is restored, "Err37: Power Failure" is displayed in the message display field.	-----	-----	<ul style="list-style-type: none"> <li>• During power failure, outer door is electrically locked.</li> <li>• After power is restored, UV resolution is repeated.</li> </ul>
	The outer door is open.	"Err38:Door Lock Failure" is displayed in the message display field.	Continuous tone (when the outer door is open)	ON (when the outer door is open)	Resolution time is extended.

**Table 9. Alarm and safety functions for STD gas calibration**

Alarm or safety function	Conditions	Display	Buzzer	Remote alarm	Safety operation
System check error at start of STD Gas calibration	The gas pressure of CO <sub>2</sub> gas line for the standard gas calibration is less than the specified pressure.	"STD Gas Calibration Error" and "Err41: STD Gas Empty" are displayed in the centre of the LCD touch panel.	Intermittent tone	ON	The STD Gas calibration is cancelled.

# TROUBLESHOOTING

If the incubator does not seem to be working properly, check the following items before calling a service engineer:

Symptom	Items to check and countermeasures
The incubator does not operate at all.	<ul style="list-style-type: none"> <li>• Is the incubator plugged in?</li> <li>• Is there a power outage, or has a circuit breaker interrupted the power?</li> <li>• The removable power supply cord is connected to the port attached on the lower right side of the cabinet.</li> </ul>
An alarm is activated.	<p>[When starting operation]</p> <ul style="list-style-type: none"> <li>• Does the chamber temperature match the set value?</li> <li>• Does the CO<sub>2</sub> gas density in the chamber match the set value?</li> </ul> <p>(1) Is the secondary pressure for the gas regulator at the specified value of 0.03 MPa(G)–0.1 MPa(G) (0.3 kgf/cm<sup>2</sup>(G)–1 kgf/cm<sup>2</sup>(G), 4.4 psi(G)–14.5 psi(G))?</p> <p>(2) Is the gas tube properly connected?</p> <p>[During operation]</p> <ul style="list-style-type: none"> <li>• Is the high-limit alarm temperature set at least 1 °C higher than the chamber set temperature?</li> <li>• Has the temperature setting been changed? Has the outer door been left open for a long time? Has a low-temperature object been placed in the chamber? If any of these is the case, the alarm will clear automatically after a short time.</li> <li>• Has the gas tube come loose, or is there a gas leak?</li> <li>• Has the CO<sub>2</sub> gas density setting been changed?</li> <li>• Is the gas cylinder empty? Check the primary pressure of the gas cylinder once a week (A primary pressure reading of 3.8 MPa(G) (38 kgf/cm<sup>2</sup>(G), 551 psi(G)) or lower is a sign that there is little gas remaining. Replace the cylinder soon).</li> <li>• Is the incubator operating beside an appliance that generates electromagnetic waves?</li> </ul>
The chamber temperature does not match the set value.	<ul style="list-style-type: none"> <li>• Is the ambient temperature less than 5 °C different from the set value for the chamber temperature?</li> <li>• Is the outer door closed with the inner door left open?</li> <li>• Is the incubator operating beside an appliance that generates electromagnetic waves?</li> </ul>
The chamber humidity does not rise.	<ul style="list-style-type: none"> <li>• Is there enough water in the humidifying pan? (Be sure to use sterile distilled water)</li> </ul>
The CO <sub>2</sub> gas density does not match the set value.	<ul style="list-style-type: none"> <li>• Is the secondary pressure for the gas regulator at the specified value of 0.03 MPa(G)–0.1 MPa(G) (0.3 kgf/cm<sup>2</sup>(G)–1 kgf/cm<sup>2</sup>(G), 4.4 psi(G)–14.5 psi(G))?</li> <li>• Is the gas tube blocked?</li> <li>• Is the duct securely attached? (See page 25)</li> <li>• Is the fan attached properly? (See page 25)</li> <li>• Is the incubator operating beside an appliance that generates electromagnetic waves?</li> </ul>
A large quantity of CO <sub>2</sub> gas is being consumed.	<ul style="list-style-type: none"> <li>• Are the outer and inner doors being frequently opened and closed?</li> <li>• Check whether gas is leaking from connectors due to deterioration of the gas tube, or whether there may be any pinhole leaks. The gas tube is a replaceable part, and it is recommended that it be replaced once a year.</li> <li>• Is the packing seal for the inner door defective?</li> <li>• Is the access hole open?</li> </ul>

Symptom	Items to check and countermeasures
Normal cultures are not possible, and the CO <sub>2</sub> gas density is suspect.	<ul style="list-style-type: none"> <li>• Is the ambient air environment around the incubator normal? Is there a source of polluted gas in the vicinity?</li> </ul>
CO <sub>2</sub> gas is not being injected.	<ul style="list-style-type: none"> <li>• CO<sub>2</sub> control for the incubator follows the ON-OFF method. CO<sub>2</sub> gas is intermittently injected as the gas density in the chamber approaches the set value. Injections may be stopped for periods of approximately 15 seconds, but that is not an error.</li> <li>• The gas is not injected until the temperature of the CO<sub>2</sub> sensor becomes sufficiently stable, approximately 1 hour after turning ON the power switch or recovering from power failure.</li> </ul>
The CO <sub>2</sub> gas density is slow to recover.	<ul style="list-style-type: none"> <li>• A HEPA filter is used for the incubator CO<sub>2</sub> gas piping. If gas density is slow to recover when the CO<sub>2</sub> gas pressure is normal, it is possible that the HEPA filter may be clogged. Contact our sales representative or agent.</li> <li>• Is there little gas remaining in the CO<sub>2</sub> gas cylinder?</li> <li>• Is the secondary pressure for the gas regulator at the specified value of 0.03 MPa(G)–0.1 MPa(G) (0.3 kgf/cm<sup>2</sup>(G)–1 kgf/cm<sup>2</sup>(G), 4.4 psi(G)–14.5 psi(G))?</li> <li>• Is the gas tube blocked?</li> <li>• Is the duct securely attached? (See page 25)</li> <li>• Is the fan attached properly? (See page 25)</li> </ul>
UV lamp lights when the outer door is open.	<ul style="list-style-type: none"> <li>• Is something pushing against the door switch?</li> </ul>
The outer door does not open.	<p>For model MCO-170AICUVH, or when the optional MCO-170EL is installed:</p> <ul style="list-style-type: none"> <li>• When the power switch is OFF, the electric lock is locked and the outer door does not open. Either turn ON the power switch or use the unlock key to override the electric lock.</li> <li>• During decontamination the outer door is electrically locked and does not open.</li> </ul>
H <sub>2</sub> O <sub>2</sub> decontamination cannot be performed.	<ul style="list-style-type: none"> <li>• If the MCO-170AICUVH is being used, is the optional MCO-HP installed?</li> <li>• If the MCO-170AICUV is being used, are the optional MCO-HP, the optional MCO-170HB, and the MCO-170EL installed?</li> <li>• If the MCO-170AIC is being used, are the optional MCO-170UVS, the MCO-HP, the optional MCO-170HB, and the MCO-170EL installed?</li> <li>• Is the UV lamp burned out?</li> <li>• Is the H<sub>2</sub>O<sub>2</sub> generator cable properly connected?</li> <li>• Has the entire bottle of the H<sub>2</sub>O<sub>2</sub> reagent been used up?</li> </ul>
In H <sub>2</sub> O <sub>2</sub> decontamination cycle, “Err34:H2O2 Volume” is activated.	<ul style="list-style-type: none"> <li>• Is the duct securely attached? (See page 25)</li> <li>• Is the fan properly attached? (See page 25)</li> <li>• Is the H<sub>2</sub>O<sub>2</sub> generator securely installed? (See page 51)</li> <li>• Has the H<sub>2</sub>O<sub>2</sub> generator reached the end of its life? If the total operating time exceeds 5,000 hours, replace the H<sub>2</sub>O<sub>2</sub> generator.</li> </ul>

# SPECIFICATIONS

Product name	CO <sub>2</sub> Incubator MCO-170AIC	CO <sub>2</sub> Incubator MCO-170AICUV	CO <sub>2</sub> Incubator MCO-170AICUVH
Medical purpose	Culture of cell tissues, organs, embryos.		
External dimensions	W 620 mm x D 730 mm x H 905 mm (W 24.4 in. x D 28.7 in. x H 35.6 in.)		
Internal dimensions	W 490 mm x D 523 mm x H 665 mm (W 19.3 in. x D 20.6 in. x H 26.2 in.)		
Interior volume	165L (5.83 cu.ft.)		
Exterior	Painted steel (Rear cover has no paint)		
Interior	Stainless steel containing copper		
Outer door	Painted steel		
Inner door	Tempered glass		
Trays	4 trays made of stainless steel containing copper W 470 mm x D 450 mm x H 12 mm (W 18.5 in. x D 17.7 in. x H 0.47 in.) Maximum load: 7 kg/tray		
Access port	Inner diameter: 30 mm (1.18 in.); on the back side		
Insulation	Styrene AcryloNitrile copolymer		
Heating system	DHA system (heater jacket + air jacket system)		
Heater	295 W		
Humidifying system	Natural evaporation with humidifying pan		
Temperature controller	PID control system		
Temperature display	Digital display		
CO <sub>2</sub> controller	PID control system		
CO <sub>2</sub> density display	Digital display		
Air circulation	Fan-assisted		
Air filter	0.3 μm, Efficiency: 99.97 % or higher		
UV lamp	-----	4 W x 1 (ozone-free emission)	
Alarms	Automatic set temperature alarm, Automatic set CO <sub>2</sub> density alarm, High-limit temperature alarm, CO <sub>2</sub> gas, various sensor/heater alarms		
Remote alarm contacts	Allowable contact capacity: DC 30 V, 2 A		
CO <sub>2</sub> inlet connection	4 mm to 6 mm (0.157 in. to 0.236 in.) diameter tube can be connected		
CO <sub>2</sub> inlet pressure	0.03 MPa(G)–0.1 MPa(G) (0.3 kgf/cm <sup>2</sup> (G)–1 kgf/cm <sup>2</sup> (G), 4.4 psi(G)–14.5 psi(G))		
Weight	80 kg		
Accessories	-----	2 unlock keys	
	1 power supply cord cover plate, 4 trays, 1 gas tube, 1 humidifying pan, 2 tube bands 1 removable power supply cord for UK 1 removable power supply cord for EU countries other than UK		

Product name	CO <sub>2</sub> Incubator MCO-170AIC	CO <sub>2</sub> Incubator MCO-170AICUV	CO <sub>2</sub> Incubator MCO-170AICUVH
Optional accessories (Refer to Table 10)	UV system set (MCO-170UVS)	standard equipment	
	H <sub>2</sub> O <sub>2</sub> decon board (MCO-170HB) Electric lock (MCO-170EL)	standard equipment	
	H <sub>2</sub> O <sub>2</sub> generator (MCO-HP)		
Optional accessories (Refer to Table 11)	Double-stacking bracket (MCO-170PS) Stacking plate (MCO-170SB)		
Optional accessories	H <sub>2</sub> O <sub>2</sub> reagent (MCO-H2O2) Gas regulator (MCO-010R) Gas auto-changer (MCO-21GC) STD gas auto-calibration kit (MCO-SG) Tray (MCO-170ST: same as that of standard accessory) Half tray (MCO-25ST) Roller base (MCO-170RB) Interface board (MCO-420MA) Small door (MCO-170ID) Interface board (MTR-L03)*; for LAN Interface board (MTR-480)*; for RS-232C/RS-485		

\* Only for the Data acquisition system MTR-5000 user.

**Note:** Refer to the updated catalogue when ordering an optional component.

Designs and specifications are subject to change without notice.

**Table 10. Optional accessories for each function**

	MCO-170AIC	MCO-170AICUV	MCO-170AICUVH
To disinfect by UV	UV system set (MCO-170UVS)	standard equipment	
To decontaminate by H <sub>2</sub> O <sub>2</sub> .	UV system set (MCO-170UVS) H <sub>2</sub> O <sub>2</sub> generator (MCO-HP) H <sub>2</sub> O <sub>2</sub> decon board (MCO-170HB) Electric lock (MCO-170EL)	H <sub>2</sub> O <sub>2</sub> generator (MCO-HP) H <sub>2</sub> O <sub>2</sub> decon board (MCO-170HB) Electric lock (MCO-170EL)	H <sub>2</sub> O <sub>2</sub> generator (MCO-HP)
To lock the outer door	Electric lock (MCO-170EL)		standard equipment

**Table 11. Bracket/plate required for each combination of double-stacking incubators**

Upper product	MCO-170AIC series MCO-170M series MCO-170AICD series			
Lower product	MCO-170AIC series MCO-170M series MCO-170AICD series	MCO-19AIC series MCO-19M series	MCO-18AC MCO-20AIC	MCO-230AIC series
Bracket/ Plate	Double-stacking bracket MCO-170PS	Stacking plate MCO-170SB		Stacking plate MCO-230SB

# PERFORMANCE

Product name	CO <sub>2</sub> Incubator MCO-170AICUVH MCO-170AICUV MCO-170AIC
Model number	MCO-170AICUVH-PE MCO-170AICUV-PE MCO-170AIC-PE
Temperature control range	Ambient temperature+5 °C to 50 °C* (ambient temperature: 5 °C to 35 °C)
Temperature distribution	±0.25 °C (ambient temperature: 23 °C, setting: 37 °C, CO <sub>2</sub> : 5 %, no load)
Temperature variation	±0.1 °C (ambient temperature: 23 °C, setting: 37 °C, CO <sub>2</sub> : 5 %, no load)
CO <sub>2</sub> control range	0 % to 20 %
CO <sub>2</sub> variation	±0.15 % (ambient temperature: 23 °C, setting: 37 °C, CO <sub>2</sub> : 5 %, no load)
Chamber humidity	95 %R.H.±5 %R.H.
Applicable environmental condition	Temperature: 5 °C to 35 °C, Humidity: 80 %R.H. max. (The optimum performance may not be obtained if the ambient temperature is not above 15 °C)
Noise level	29 dB (A scale)
Power consumption	Max. 379 W
Heat emission	Max. 1,070 kJ/h
Rated voltage, frequency	AC 220 V–240 V, 50 Hz/60 Hz
Amperage	Max. 1.8 A

\* When set temperature is 37 °C, ambient temperature must be 32 °C or less. Regardless of ambient temperature, the maximum of temperature control range is always 50 °C.

**Note:** The unit with CE mark complies with EU directives.  
Based on our measuring method.

## ⚠ CAUTION

**Please fill in this form before servicing.  
Hand this form to the service engineer to keep for his and your safety.**

### Safety check sheet

1. Unit contents :

Risk of infection: Yes No  
 Risk of toxicity: Yes No  
 Risk from radioactive sources: Yes No

(List all potentially hazardous materials that have been stored in this unit.)

Notes :

2. Contamination of the unit

Unit interior

No contamination Yes No

Decontaminated Yes No

Contaminated Yes No

Others:

3. Instructions for safe repair/maintenance/disposal of the unit

a) The unit is safe Yes No

b) There is some danger (see below) Yes No

Procedure to be adhered to in order to reduce safety risk indicated in b) below.

Date :

Signature :

Address, Division :

Telephone :

Product name : CO <sub>2</sub> Incubator	Model No. MCO-	Serial number :	Date of Installation :
---	-------------------	-----------------	------------------------

Please decontaminate the unit yourself before calling the service engineer.

# DISPOSAL OF UNIT

Before disposal of this CO<sub>2</sub> incubator, contact our sales representative or agent for further information. Improper handling of biohazardous waste can result in accidental exposure to infectious agents. If there is a danger of biohazard, decontaminate the incubator as thoroughly as possible before disposal.

## Disposal of Old Equipment and Batteries Only for European Union and countries with recycling systems

These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries must not be mixed with general household waste.

For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points in accordance with your national legislation.

By disposing of them correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment.

For more information about collection and recycling, please contact your local municipality.

Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

### Note for the battery symbol (bottom symbol):


This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

Label indication is obliged to comply with Taiwanese battery regulation.

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The screenshot shows the 'Key Lock' configuration screen. At the top, there is a title 'Key Lock' with a square icon to its left. To the right of the title are two buttons: '↑ Top' and '← Back'. Below the title, there are three rows of controls:

- The first row is 'Key Lock :', followed by a slider control currently set to 'ON'.
- The second row is 'Password :', followed by a text input field containing four asterisks '\*\*\*\*'.
- The third row is 'Confirm Password :', followed by another text input field containing four asterisks '\*\*\*\*'.

At the bottom right of the screen is an 'Apply' button.

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The screenshot shows the 'Key Lock' configuration screen. At the top, there is a title 'Key Lock' with a square icon to its left. To the right of the title are two buttons: '↑ Top' and '← Back'. Below the title, there are three rows of controls:

- The first row is 'Key Lock :', followed by a slider control currently set to 'OFF'.
- The second row is 'Password :', followed by an empty text input field.
- The third row is 'Confirm Password :', followed by another empty text input field.

At the bottom right of the screen is an 'Apply' button.

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

The screenshot shows the 'Key Lock' configuration screen. At the top right, there are '↑ Top' and '← Back' buttons. The 'Key Lock' status is set to 'ON'. Below this, there are two password input fields: 'Password' and 'Confirm Password', both containing '\*\*\*\*'. An 'Apply' button is located at the bottom right.

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The screenshot shows the 'Key Lock' configuration screen. At the top right, there are '↑ Top' and '← Back' buttons. The 'Key Lock' status is set to 'OFF'. Below this, there are two empty password input fields: 'Password' and 'Confirm Password'. An 'Apply' button is located at the bottom right.

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

The screenshot shows a dark-themed menu titled "Key Lock". At the top right are two buttons: "↑ Top" and "← Back". Below the title, there are three input fields. The first is labeled "Key Lock" and contains a slider control with a white bar on the left and the text "OFF" on the right. The second is labeled "Password" and contains a grey rectangular field. The third is labeled "Confirm Password:" and contains a grey rectangular field. At the bottom right is a button labeled "Apply".

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The screenshot shows the same "Key Lock" menu. The slider control now has a white bar on the right and the text "ON" on the left. The "Password" field contains four asterisks "\*\*\*\*". The "Confirm Password:" field also contains four asterisks "\*\*\*\*". The "Apply" button remains at the bottom right.

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The screenshot shows a dark-themed control interface titled "Stand-by Setting". At the top right are two buttons: "↑ Top" and "◀ Back". The main area contains three rows of settings:

- Temperature : 37.0 °C
- CO<sub>2</sub> Density : 5.0 % with a toggle switch set to "ON"
- High Limit : 53.0 °C

An "Apply" button is located at the bottom right of the screen.

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

The screenshot shows the 'Key Lock' configuration screen. At the top, there is a title 'Key Lock' with a square icon to its left. In the top right corner, there are two buttons: '↑ Top' and '◀ Back'. Below the title, there are three input fields. The first is 'Key Lock' with a toggle switch currently set to 'ON'. The second is 'Password' with a text box containing five asterisks. The third is 'Confirm Password' with a text box also containing five asterisks. At the bottom right, there is an 'Apply' button.

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The screenshot shows the 'Key Lock' configuration screen. At the top, there is a title 'Key Lock' with a square icon to its left. In the top right corner, there are two buttons: '↑ Top' and '◀ Back'. Below the title, there are three input fields. The first is 'Key Lock' with a toggle switch currently set to 'OFF'. The second is 'Password' with an empty text box. The third is 'Confirm Password' with an empty text box. At the bottom right, there is an 'Apply' button.

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

The screenshot shows a dark-themed menu titled "Key Lock". At the top right, there are two buttons: "↑ Top" and "◀ Back". The "Key Lock" label is followed by a colon and a slider control. The slider is positioned to the left, and the word "OFF" is displayed in a white box on the slider. Below this, there are two input fields for "Password" and "Confirm Password", both containing greyed-out text. At the bottom right, there is an "Apply" button.

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
The screenshot shows the same "Key Lock" menu. The "Key Lock" label is followed by a colon and a slider control. The slider is positioned to the right, and the word "ON" is displayed in a white box on the slider. Below this, the "Password" and "Confirm Password" input fields now contain four asterisks "\*\*\*\*". The "Apply" button remains at the bottom right.

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

The screenshot shows the 'Key Lock' configuration screen. At the top left is a square icon and the text 'Key Lock'. At the top right are two buttons: '↑ Top' and '← Back'. Below this, there are three rows of input fields. The first row is 'Key Lock :', followed by a toggle switch currently set to 'ON'. The second row is 'Password :', followed by a text input field containing six asterisks. The third row is 'Confirm Password :', followed by another text input field containing six asterisks. At the bottom right is an 'Apply' button.

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The screenshot shows the 'Key Lock' configuration screen. At the top left is a square icon and the text 'Key Lock'. At the top right are two buttons: '↑ Top' and '← Back'. Below this, there are three rows of input fields. The first row is 'Key Lock :', followed by a toggle switch currently set to 'OFF'. The second row is 'Password :', followed by a text input field that is currently empty. The third row is 'Confirm Password :', followed by another text input field that is currently empty. At the bottom right is an 'Apply' button.

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**PHC Europe B.V.**