EVERmed

MEDICAL REFRIGERATION 

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User Instruction

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ARMADIO REFRIGERATO

REFRIGERATED CABINET

ARMOIRE REFRIGEREE

KÜHLSCHRANK

STANDARD

VERS. 0816

CE

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ENGLISH
Content
Section 8 USER INSTRUCTIONS

8.1 CONTROLS	
8.1.1 Description of controls and pushbuttons.....	10
8.1.1.1 Keyboard safety locking protection	10
8.1.2 Visual functioning indications	10
8.2 OPERATION OF THE CONTROL PANEL.....	11
8.2.1 Start-up	11
8.2.2 Stopping the unit	11
8.2.3 Set-up and adjustment operations	11
8.2.4 Automatic and manual defrosting	11
8.2.5 Alarm signallings	11
8.2.6 Alarm memorization.....	12
8.2.7 Display the alarm memorized.....	12
8.2.8 Reset the alarm list.....	12
8.2.9 Compressor working hours	12
8.2.10 Display the compressor working hours.....	12
8.2.11 Reset the compressor working hours	12
8.3 OPERATIONS OF THE TEMPERATURE CHART RECORDER.....	13
8.3.1 Description of the temperature chart recorder	13
8.3.1.1 Paper chart replacement	13
8.3.1.2 Ink-tip replacement	13
8.3.1.3 Battery replacement	13
8.4 OPERATION OF THE DIGITAL PRINTER	14
8.4.1 Description of the digital printer	14

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Section 8 USER INSTRUCTIONS

The information in this section of the manual regards the user or other non-specialized personnel (see par. 1.3 in the “Instruction and Maintenance Manual”). After the appliance has been installed in accordance with the instructions of section 3 of this manual, it is ready for use.


8.1 CONTROLS

According to the models, the appliance is equipped of different types of controls:

- Control panel (Fig. 40)
- Temperature chart recorder (Fig. 11)
- Digital printer

8.1.1 Description of controls and pushbuttons (Fig.40)

Fig.40

The internal temperature is automatically displayed and controlled through a digital control thermo-regulator and the lighting icons give indication of the functioning status. The temperature controller is provided with 6 keys with the following specific functions:

- Display (11), to check temperature and operation of the machine.
 - STAND-BY pushbutton (1). If pressed for 5 seconds light on the thermo-regulator
 - SET pushbutton (13). If pressed allows the regulation of the working set of the machine
 - UP pushbutton (6). It allows the increasing of the values, (higher temperature or higher values during the programming phase). If pressed, signs on the display an alarm indicating “rLS”.
 - DOWN pushbutton (5). If pressed allows the decrease of the values (lower temperature or lower values during programming phase). If pressed, for more than 2 seconds, accede directly to the under menu of probes temperature, alarm visualization, compressor working hours and regulation of the internal watch.
 - DEFROSTING pushbutton (16) progress a manual defrosting cycle.
 - LIGHTING pushbutton (17) lights the internal lighting .
- All the keys, if pressed, have also the function of silencing the buzzer of the thermo-regulator and memorize the possible alarms.

8.1.1.1 Keyboard safety locking protection

The control panel keyboard is provided with a keys locking protection, to avoid tampering from unauthorized personnel. The key locking protection has to be activated from the user if necessary.

- press at the same time DOWN and ON/STANDBY pushbuttons for 2 seconds and the display will show “Loc” (keyboard locked). If the keyboard is locked, it will not be possible to execute any operation with the keys and if pressed the display will show “Loc”.
- To unlock the keyboard:
- press at the same time DOWN and ON/STANDBY pushbuttons for 2 seconds and the display will show “UnL” (keyboard unlocked)
- At this stage it will be possible to use again the keys and their functionalities.

8.1.2 Visual functioning indications (Fig. 40)

The thermo-regulator lights the following icons for functions in progress:

- ON/STANDBY Led (20), turns on when the machine is connected but in stand by.
- COMP led (10) lights steadily when the compressor is running, blinks in case of starting delay or when protection is active
- FAN led (8) lights steadily when the evaporator fan is working
- DEF led (9) lights steadily during defrosting cycle
- ALL led (7) lights steadily for alarm of temperature and probe faulty
- HACCP led (14) lights steadily or intermittent to signalized alarm of an HACCP memorized

In the thermo-regulator appears the followings signals for functions in progress:

- CH led (15) (adjustable wrench) lightly or intermittent signalizes a request of program maintenance for the compressor working hours.
- LI led (19) lightly when the internal light is on

8.2 OPERATION OF THE STANDARD CONTROL PANEL

8.2.1 Start-up

Before starting up the unit check that the electrical connections have been made correctly as indicated in headings 3.3 and 3.4 "Instruction and Maintenance Manual".

Perform preliminary cleaning of the unit as described in sub-heading 5.2.1. "Instruction and Maintenance Manual"

Starting sequence (Fig. 40)

- press STAND-BY pushbutton for 5 seconds and the display lights up
- wait until the display stop blinking
- set the internal watch with the actual hour

The first step to do is the regulation of the internal watch in the following way:

- press the DOWN key for 2 seconds. The display will evidence "rtC" (real time clock)
- press the SET key and the display will evidence yy00 (year). Press up or down to change it
- press the SET key and the display will show up nn01 (month). Press up or down to change it
- press the SET key and the display will show up dd01 (day). Press up or down to change it
- press the SET key and the display will show up hh00 (hour). Press up or down to change it
- press the SET key and the display will show up nn00 (minutes). Press up or down to change it and then press SET.

At the end of the complete set up do not operate on the pushbuttons for 60 seconds.

Automatically it will memorize the new settings and exit of the process.

8.2.2 Stopping the unit (Fig.40)

Press STAND-BY pushbutton for 5 seconds and the display switches off. In case the appliance will remain off for a prolonged period, it is strongly recommendable to remove the appliance plug from the power socket.

8.2.3 Set-up and adjustment operations

The appliance is factory set for operation at the following temperature:

- serie MPR : +2°C / +15°C
- serie LR : 0°C / +15°C
- serie LF : -5°C / -20°C ~ -25°C
- serie LDF : -15°C / -30°C
- serie PDF : -25°C / -40°C
- serie MPRR : +2°C / +15°C (1st refrigerator compartment), +2°C / +15°C (2nd refrigerator compartment)
- serie LCRR : 0°C / +15°C (1st refrigerator compartment), 0°C / +15°C (2nd refrigerator compartment)
- serie LCRF : 0°C / +15°C (refrigerator compartment), -5°C / -20°C ~ -25°C (freezer compartment)

If you intend to use different temperatures, proceed as follows (Fig.40):

- press SET (13) and release. The display shows the actual programmed temperature
- to increase the temperature, press the UP key (6) within 5 seconds
- to lower the temperature, press the DOWN key (5) within 5 seconds

The instrument automatically memorizes the programmed temperature last value.

8.2.4 Automatic and manual defrosting

The appliance is factory set for automatic defrosting at pre-set intervals.

If you wish to defrost manually in accordance with effective requirements, proceed as follows (Fig. 40):

- press the DEFROST pushbutton (16) for more than 3 seconds

During the automatic or manual defrost cycle the DEF led (9) on the control panel will be steadily illuminated. At the end of the defrost cycle the led will switch off and the appliance automatically resumes its normal operating cycle.

8.2.5 Alarm signalings

The thermo-regulator, in occasion of an alarm condition, will activate a buzzer for the acoustic signalling and light-up the led 7, and it will signal, with a code indication, the type of alarm. The display will indicate:

Alarm code	Description
Pr1	Faulty chamber probe
Pr2	Faulty evaporator probe
Pr3	Faulty condenser probe
COH	Overheated condenser
DFd	End of defrosting with after maximum time allowed
AL	Low temperature alarm
AH	High temperature alarm
PF	Power failure
id	Door ajar alarm
Rtc	Wrong programmation of the internal watch (must be programmed again)
RLS	Signals the presence of an alarm in the memory
LS	Alarm folder

The thermo-regulator can store, into its memory, up to 3 alarms (HACCP). The controller provide the following information :

- critical value of temperature
- date and real time of alarm registration
- duration of the alarm (from 1 until 99h e 59 min., partial if alarm in progress)

The values of the alarm limits are factory set and change according to the product series.

8.2.6 Alarm memorization

When the value of the temperature of the appliance goes out of the pre-fixed minimum or maximum limits, it will appear an alarm signal and automatically create a folder "LS" in the thermo-regulator menu. The folder contains the highest and lowest values achieved from the temperature and the time of the alarm in progress registered.

8.2.7 Display the alarm memorized

The alarms generated by the thermo-regulator are displayed as follows:

- press DOWN key (5), for 2 seconds. The display will show the first label available
- press UP key (6) or DOWN (5) to select the "LS" label (alarms folder)
- press SET key (13). The display will evidence the type of alarm (AL, AH)
- press once more the SET key (13) to see the temperature value and alarm duration, moreover date and real time of the alarm registration and conditions in progress.

Example:

AH1	High temperature alarm
20	Highest temperature achieved
Sta	Viewer of date and hour of the alarm
y (11)	Year of alarm registration
n (09)	Month of alarm registration
d (15)	Day of alarm registration
h (16)	Hour of alarm registration
n (30)	Minutes of alarm registration
Dur	Viewer of the alarm duration
h (2)	Hour of duration of the alarm conditions
n (30)	Minutes of duration of the alarm conditions

In the above example, the thermo-regulator registered an alarm for high temperature (AH1), the temperature has reached 20°C on the 15th of September 2011 at 04.30 pm and the alarm lasted 2 h and 30 minutes.

To the exit of the alarm displaying don't operate in any other key for 15 seconds.

8.2.8 Reset of the alarm list

To reset the alarm folders proceed in the following way:

- press DOWN pushbutton (5) for 2 seconds
- press UP or DOWN pushbutton to select label "rLS"
- press SET pushbutton (13), and see 0
- press UP pushbutton (6) within 15 minutes to set the value 149
- press SET pushbutton (13) and don't act any other operation kind for 15 seconds

The display will show ----- in intermittance for 4 seconds and the icon HACCP will turn off. If there isn't any other kind of alarm in memory, the label "rLS" will not be shown up. If the alarm folder will not be delete the most recent alarm will substitute the oldest one.

8.2.9 Compressor working hours

The controller is able to memorize until 9999 of the compressor working hours to program the maintenance operation.

At the expiry of the programmed hours suggested as maintenance interval, the icon 15 (wrench) will light up on the display.

8.2.10 Display the compressor working hours

- press for 2 seconds the DOWN pushbutton (5) and the display will show the first label available
- press UP or DOWN pushbuttons to select CH
- press SET and the display will show the working hours of the compressor

To exit press SET (13) or don't operate for the next 15 seconds.

8.2.11 Reset the compressor working hours

- press for 2 seconds the DOWN pushbutton (5) and the display will show the first label available
- press UP or DOWN pushbuttons to select "rCH"
- press SET (13) and the display will show 0 (zero)
- press UP pushbutton (6) and program 149
- press SET (13) and don't operate for the next 15 seconds. The display will show ----- blinking for 4 seconds and the value of the hours will be deleted.

For more instructions regarding the thermo-regulator function see the manual in attachment.

8.3 OPERATION OF THE TEMPERATURE CHART RECORDER

8.3.1 Description of the temperature chart recorder

The appliance may be fitted (present if ordered in some series), with the temperature chart recorder, recording on paper charts the appliance internal temperature. The temperature recorder is available in different temperature ranges according to the model on which it is installed:

- Range : -10°C / +40°C
- Range : -35°C / +15°C
- Range : -100°C / +50°C

The survey of the chamber temperature is done with a sensor independent from the control panel sensor, while the recording, weekly operation, is on paper charts with an ink-tip. The temperature recording is granted even during power failure periods thanks to the battery powering of the recording device. The battery is 1,5V, AA type. A perspex cover allows checking the diagrams but at the same time it protects them against tampering of unauthorized personnel thanks to the key locks safety.

8.3.1.1 Paper chart replacement

- Open the perspex cover
- Unscrew the fixing nut
- Lift gently the metal arm on which is fitted the ink-tip
- Remove the installed paper chart, paying attention at the 2 metal edges surrounding the chart that keep it in position
- Insert the new paper chart, paying attention to insert it into the central hinge and into the 2 metal edges surrounding the chart
- Lower gently the metal arm on which is fitted the ink-tip in order to return it in the original position
- Drive the paper chart in order to position the ink-tip writing point in the exact point from which it has to start to record the temperature, paying attention to do not write on the diagram. In order to choose the exact recording start point, refer to the days and times printed on the chart itself
- Screw the fixing nut
- Close the perspex cover





8.3.1.2 Ink-tip replacement

- Open the perspex cover
- Lift gently the metal arm on which is fitted the ink-tip
- Remove, pulling, the ink-tip from the metal arm
- Insert the new ink-tip on the metal arm paying attention to the guideways on the ink-tip itself. Push until reach the limit-stop
- Lower gently the metal arm on which is fitted the ink-tip in order to return it in the original position
- Close the perspex cover

8.3.1.3 Battery replacement

- Open the perspex cover
- Unscrew the fixing nut
- Lift gently the metal arm on which is fitted the ink-tip
- Remove the installed paper chart, paying attention at the 2 metal edges surrounding the chart that keep it in position
- Pull the central hinge on which was screwed the nut, and at the same time moving slightly up and down to make easier the extraction of the clockwise mechanism
- Replace the battery in the rear side of the clockwise mechanism paying attention to the battery polarity
- Insert the clockwise mechanism in the hole left on the chart recorder, simply making pressure on the mechanism
- Insert the paper chart, paying attention to insert it into the central hinge and into the 2 metal edges surrounding the chart
- Lower gently the metal arm on which is fitted the ink-tip in order to return it in the original position
- Drive the paper chart in order to position the ink-tip writing point in the exact point from which it has to start to record the temperature, paying attention to do not write on the diagram. In order to choose the exact recording start point, refer to the days and times printed on the chart itself
- Screw the fixing nut
- Close the perspex cover

Fig.11

8.4 OPERATION OF THE DIGITAL PRINTER

8.4.1 Description of the digital printer

The appliance may be fitted (present if ordered in some series), with the temperature digital printer.

The digital printer records the temperature in an internal non-volatile memory the temperature of the appliance. The recording periods and intervals are fully adjustable and can be easily set through the printer soft-touchpads and the wide display (with clear written indications). Then, the temperature trends can be printed on paper, selecting among 3 printing modes according to the need, online (printing out at fix intervals), on daily basis (once per day) or an historical report of all the datas acquired. The survey of the chamber temperature is done with a sensor independent from the control panel sensor. To increase the safety level, it is included an alarm for sensor failure. The printing is ink-free on thermal paper.

For more detailed information, please refer to the specific booklet (included among the instructions if the printer is present).

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EVERmed Srl: Via Galileo Galilei, 2 - 46020 Motteggiana (MN) ITALY